

Front Riggings Index and Calf Pads

Due to a change to the Legrest Calf Pad size made on 4/23/12, you must also order the hardware kit (1066591) when ordering replacement Calf Pads.

Calf Pads

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	B,H	88804530X003	Pad, Legrest, Vinyl - Right	1	1
	B,H	88804530X004	Pad, Legrest, Vinyl - Left	1	1
	A,C,H	88804529X003	Pad, Legrest, Vinyl - Right	1	1
	A,C,H	88804529X004	Pad, Legrest, Vinyl - Left	1	1
	A,D,E, H	1091767-NLA	NLA - Pad with Screws, Legrest (Pair)	2	2
	A,D,F	8881091759	Pad w/ Cover and Screws, Legrest Vinyl (Pair)	2	2
	A,D,F, H	8881091764	Pad, Legrest, Vinyl - Right/Left	1	2
	G	1066591	Hardware, Legrest Pad	1	1
			Screw, Phillips Pan Head (#10-32 x 11/16")	1	2
			Screw, Phillips Pan Head (#10-32 x 1")	1	2
			Instruction Sheet, Legrest Pad Mtg.	1	1
	A,D	01415X000-NLA	NLA - Cover, Legrest Pad	1	1
		See Chart	Front Riggings Index	1	1
NOTE: A - Specify upholstery color, see page titled "Frames Finishes and Upholstery Colors" B - 12" Calf Pads C - 14-15" Calf Pads D - 16-22" E - Before 9/5/00 F - After 9/4/00 G - Includes Screws and Instruction Sheet. H - Due to a change to the Legrest Calf Pad size made on 4/23/12, you must also order the hardware kit (1066591) when ordering replacement Calf Pads.					

Front Riggings Index

Front Rigging	Chair Width	Support Assembly Right	Support Assembly Left	Pivot Slide Tube w/Aluminum Footplate	Pivot Slide Tube w/Aluminum Footplate
				See Note(s) A	See Note(s) A
H93E	10-12"	1040715-NLA	1040716-NLA	1007274-NLA	1007275-NLA
	14"	1040715-NLA	1040716-NLA	1033506-NLA	1033507-NLA
H93M	10-12"	1035172-NLA	1035171-NLA	1007274-NLA	1007275-NLA
	14"	1035172-NLA	1035171-NLA	1033506-NLA	1033507-NLA
H903MC	12"	1031194-NLA	1031195-NLA	1007274-NLA	1007275-NLA
	14"	1031194-NLA	1031195-NLA	1033506-NLA	1033507-NLA
H903MK	12"	8881031198-NLA	8881031199-NLA	1007274-NLA	1007275-NLA
	14"	8881031198-NLA	8881031199-NLA	1033506-NLA	1033507-NLA
H990C	10"	1020011-NLA See Note(s) B	1020012-NLA See Note(s) B	N/A	N/A
	12-14"	1020015-NLA See Note(s) B	1020016-NLA See Note(s) B	N/A	N/A
H990K	10"	1020009-NLA See Note(s) B	1020010-NLA See Note(s) B	N/A	N/A

Front Riggings Index and Calf Pads

Due to a change to the Legrest Calf Pad size made on 4/23/12, you must also order the hardware kit (1066591) when ordering replacement Calf Pads.

Front Riggings Index

Front Rigging	Chair Width	Support Assembly Right	Support Assembly Left	Pivot Slide Tube w/Aluminum Footplate See Note(s) A	Pivot Slide Tube w/Aluminum Footplate See Note(s) A
	12-14"	1020013-NLA See Note(s) B	1020014-NLA See Note(s) B	N/A	N/A
H904M	12"	8881031465-NLA	8881031466-NLA	1031471-NLA	1031472-NLA
	14"	1031467-NLA	1031468-NLA	1020234-NLA	1020235-NLA
NOTE: A - Aluminum footplates on chairs built before 2/96 had a "Ball Burnish" finish. Footplates after 1/96 now have a black finish. When replacing only 1 pivot slide tube w/footplate before 2/96, it is suggested that a second new (black) footplate be ordered. B - Pivot slide tube w/footplate is included with Support Assembly.					
Orders for Front Riggings in pairs must be placed through the Configurator. Please contact our Customer Service Coordinator.					