

Front Riggings - Calf Pads

Due to a change to the Legrest Calf Pad size made on 4/23/12, you must also order the hardware kit (1066591) when ordering replacement Calf Pads.

Calf Pads

Part Number	Description	Chair Width
88804529X003 See Note(s) C, A	Pad, Legrest, Vinyl - Right	14-15" Wide - Right Side
88804529X004 See Note(s) C, A	Pad, Legrest, Vinyl - Left	14-15" Wide - Left Side
8881091819 See Note(s) C, B	Kit, Calf Pad w/ Screws, Expanded Vinyl - Right/Left (16 - 22" W)	16-22" Wide
1091800 See Note(s) C, B	Kit, Calf Pad w/ Screws, Black Plastic - Right/Left	16-22" Wide
NOTE: A - Hardware is not included. B - Hardware included. C - Due to a change to the Legrest Calf Pad size made on 4/23/12, you must also order the hardware kit (1066591) when ordering replacement Calf Pads.		

Calf Pad Hardware

Part Number	Description
1066591 See Note(s) A	Hardware, Calf Pad
NOTE: A - Used for padded calf pads only, includes Screws and Instruction Sheet.	