




# Disinfectant System Components




1




2




3




4




5



6



7




8


10




9




11



12



13



14

## Disinfectant System Components

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	PAR1110C	Valve, Two Way Ball (Disinfectant on/off)	1	1
2		PAR1110	Handle, Replacement	1	1
3		PAR1108	Metering Valve	1	1
4		IHPAR1106	Metering Valve Replacement Knob	1	1
5		PAR5673	Check Valve (1/4" FPT)	1	1
6		PAR1030	Injector Assembly (3/8" FPT)	1	1
7		IHPAR1034	Check Valve and Suction Tube, Disinfectant	1	1
8	B	IHPAR1032-NLA	NLA - Weight, Ceramic (for disinfectant suction tube and check valve)	1	1
9		PAR7380AG	Back Flow Preventer - Conbraco (1/2") (One Screw in Cap)	1	1
10		IHPAR7383C	Bonnet	1	1
11		IHPAR7384C	Kit, Repair	1	1
12		PAR1136	Elbow, 90 Degrees - Chrome	1	1
13		PAR6160D	Washer, Nylon, Backflow Preventer and Chrome Elbow	1	1
14		IHPAR3617A-NLA	NLA - Washer, Disinfectant Burp Fitting	1	1
		IHPAR3617-NLA	NLA - Burp Fitting, Disinfectant	1	1

NOTE: A - Replaces IHPAR1111C (No Longer Available) on all Tubs built after June 2009. Replacements for Tubs built prior to June 2009 must have rinse line capped when replacing fitting.  
B - Not shown in picture.