


Seat Frame Attaching Hardware

Adult 16" - 24" Wide



Seat Frame Attaching Hardware Adult 16" - 24" Wide

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1043627-NLA	NLA - Screw, Hex Head w/ Patch (5/16-18 x 1-3/4")	1	4
1	B	1070398-NLA	NLA - Screw, Hex Head (5/16-18 x 1-1/2")	1	4
2	A	1070413-NLA	NLA - Spacer (21/64 x 1/2 x 17/64")	1	4
3		1043627-NLA	NLA - Screw, Hex Head w/ Patch (5/16-18 x 1-3/4")	1	4
4		1025704	Washer, Flat (5/16 x 5/8 x 1/16")	1	12
5		1070388	Support, Seat, Black	1	2
6		1069724	T-Nut (5/16-18)	1	5
7		1070462-NLA	NLA - Screw, Hex Head w/ Patch (5/16-18 x 7/8")	1	4
8		1026158	Locknut (5/16-18)	1	4
9		1070396-NLA	NLA - Strap, Bolt, Black	1	2
10	D	1070343-NLA	NLA - Plate, Seat Mounting Low, Right	1	2
10	D	1070344-NLA	NLA - Plate, Seat Mounting Low, Left	1	2
10	E	1070345-NLA	NLA - Plate, Seat Mounting Medium, Right	1	2
10	E	1070346-NLA	NLA - Plate, Seat Mounting Medium, Left	1	2
10	F	1070347-NLA	NLA - Plate, Seat Mounting Tall, Right	1	2
10	F	1070348-NLA	NLA - Plate, Seat Mounting Tall, Left	1	2
11	C	1093314-NLA	NLA - Hardware, Attaching, Seat Frame (16" W)	1	1
11	C	1093311-NLA	NLA - Hardware, Attaching, Seat Frame (17" and Wider)	1	1

NOTE: A - Part # 1043627 and 1070413 are used on the 16" wide chairs only.
 B - Part # 1070398 is used on 17" and wider chairs only.
 C - Item number 10 not included in seat frame attaching hardware.
 D - 18-1/4" Seat to Floor
 E - 20-1/2" Seat to Floor
 F - 21-3/4" Seat to Floor