


Crossbrace - Conventional Fixed Seat Depth



Crossbrace - Conventional Fixed Seat Depth

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	#1 Crossbrace Assembly - Junior Fixed	1	1
2	B	1173765	Kit, Saddle Clamp w/ Hardware (Pair)	1	1
3			Clamp, Saddle	1	2
4			NLA - Screw, Socket Button Head (M6 x 40mm)	1	4
5			NLA - Locknut (M6 x 1)	1	4
6	C	1173776	Kit, Pivot Link w/ Hardware (12 - 15" W)	1	1
7			Link, Pivot, Crossbrace (12-15" W)	1	2
8			Screw (M6 x 45mm)	1	2
9			NLA - Locknut (M6 x 1)	1	4
10	D	1173777	Kit, Crossbrace w/ Hardware (Pair)	1	1
11			Washer, Plastic	1	2
12			Bolt (M6 x 80mm)	1	2
13			Washer (6.6 x 22 x 2)	1	4
14			NLA - Locknut (M6 x 1)	1	2
15	E	1173774-NLA	NLA - Kit, Plug Buttons w/ Hardware	1	1
16			NLA - Plug Button	1	4
17			NLA - Screw, Phillips Pan Head (M3.5 x .6 x 9.5mm)	1	4

NOTE: A - Includes items 2 - 17.
 B - Includes items 3 - 5.
 C - Includes items 7 - 9.
 D - Includes items 11 - 14.
 E - Includes items 16 - 17.

#1 Crossbrace Assembly - Junior Fixed

Seat Depth	12" Wide Chair	13" Wide Chair	14" Wide Chair
12" Deep	1172267	N/A	N/A
13" Deep	1172268	N/A	N/A
14" Deep	1172269	1172270	1172272
15" Deep	N/A	1172271	1172273