


Seat Frame Attaching Hardware

Junior 12" - 16" Wide



Seat Frame Attaching Hardware Junior 12" - 16" Wide

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1043627-NLA	NLA - Screw, Hex Head w/ Patch (5/16-18 x 1-3/4")	1	8
2		1026158	Locknut (5/16-18)	1	8
3		1070465-NLA	NLA - Bracket, Mounting, Black - Junior	1	2
4		1070462-NLA	NLA - Screw, Hex Head w/ Patch (5/16-18 x 7/8")	1	8
5		1025704	Washer, Flat (5/16 x 5/8 x 1/16")	1	12
6		1070388	Support, Seat, Black	1	2
7		1069724	T-Nut (5/16-18)	1	5
8	B	1070343-NLA	NLA - Plate, Seat Mounting Low, Right	1	2
8	B	1070344-NLA	NLA - Plate, Seat Mounting Low, Left	1	2
8	C	1070345-NLA	NLA - Plate, Seat Mounting Medium, Right	1	2
8	C	1070346-NLA	NLA - Plate, Seat Mounting Medium, Left	1	2
8	D	1070347-NLA	NLA - Plate, Seat Mounting Tall, Right	1	2
8	D	1070348-NLA	NLA - Plate, Seat Mounting Tall, Left	1	2
9		1070396-NLA	NLA - Strap, Bolt, Black	1	2
10	A	1093313-NLA	NLA - Hardware, Seat Mounting, Junior (12-16" Wide)	1	1

NOTE: A - Item number 8 not included in seat frame attaching hardware.
 B - 18-1/4" Seat to Floor
 C - 20-1/2" Seat to Floor
 D - 21-3/4" Seat to Floor