



# Backrests and Hardware

## w/o Side Guards or Other Side Guards



## Backrests and Hardware w/o Side Guards or Other Side Guards

| Item Number                      | Note      | Part Number                | Description                                 | Quantity Per |          |
|----------------------------------|-----------|----------------------------|---------------------------------------------|--------------|----------|
|                                  |           |                            |                                             | Package      | Assembly |
| 1                                | A         | 1114096-NLA                | NLA - Kit, Fold Down Back Mounting Hardware | 1            | 1        |
| 2                                |           | 1114084-NLA                | NLA - Bracket, Support Back, Black -Right   | 1            | 1        |
| 2                                |           | 1114085-NLA                | NLA - Bracket, Support Back, Black -Left    | 1            | 1        |
| 3                                |           | 1116005-NLA                | NLA - Bracket, Latch, Black                 | 1            | 2        |
| 4                                |           | 1039232-NLA                | NLA - Cam, Adjustment                       | 1            | 2        |
| 5                                |           | 40157X020                  | Spacer, Coved (1/4 x 1 x 9/16")             | 1            | 4        |
| 6                                |           | 1040776-NLA                | NLA - Spacer, Nylon (1/4 x 1 x 1/4")        | 1            | 4        |
| 7                                |           | 1114099                    | NLA - Screw, Hex Head (1/4-28 x 2")         | 1            | 8        |
| 8                                |           | 1071496-NLA                | NLA - Pin, Roll (3/16 x 3/4")               | 1            | 2        |
| 9                                |           | 1025728                    | Washer (1/4 x 1/2 x 1/32")                  | 1            | 2        |
| 10                               |           | 1114077-NLA                | NLA - Latch, Torsion Spring - Right         | 1            | 1        |
| 10                               |           | 1114078-NLA                | NLA - Latch, Torsion Spring - Left          | 1            | 1        |
| 11                               |           | 1026162                    | Spacer (1/4 x 3/8 x 11/32")                 | 1            | 4        |
| 12                               |           | 1058636-NLA                | NLA - Washer, Nylon (25/64 x 7/8 x 1/32")   | 1            | 2        |
| 13                               |           | 1026907                    | Washer (1/4 x 1/2 x 1/16")                  | 1            | 14       |
| 14                               |           | 1025765                    | Locknut (1/4-28)                            | 1            | 10       |
| 15                               |           | 1141633                    | Cord, Nylon (6-1/2 ft.)                     | 1            | 1        |
| 16                               |           | 1039233                    | Spacer (1/4 x 7/16 x 5/8")                  | 1            | 2        |
| 17                               |           | 1045495-NLA                | NLA - Screw, Set Socket (1/4-20 x 1/4")     | 1            | 2        |
| 18                               |           | See Chart                  | #18 - One Piece Backrest Frame              | 1            | 1        |
| 19                               | See Chart | #19 - Upright Tube         | 1                                           | 2            |          |
| 20                               | See Chart | #20 - Backrest Crossmember | 1                                           | 1            |          |
| NOTE: A - Includes items 2 - 16. |           |                            |                                             |              |          |

### #18 - One Piece Backrest Frame

| Chair Width | 8-11" High  | 10-14" High | 12-16" High | 14-18" High | 16-20" High |
|-------------|-------------|-------------|-------------|-------------|-------------|
| 12"         | 1071788-NLA | 1071789-NLA | 1071790-NLA | 1071791-NLA | 1071792-NLA |
| 13"         | 1071793-NLA | 1071794-NLA | 1071795-NLA | 1071796-NLA | 1071797-NLA |
| 14"         | 1024706-NLA | 1026188-NLA | 1042982-NLA | 1026757-NLA | 1026764-NLA |
| 15"         | 1024707-NLA | 1026189-NLA | 1042986-NLA | 1026758-NLA | 1026765-NLA |
| 16"         | 1024708-NLA | 1026190-NLA | 1042987-NLA | 1026759-NLA | 1026766-NLA |
| 17"         | 1024709-NLA | 1026191-NLA | 1042988-NLA | 1026760-NLA | 1026767-NLA |
| 18"         | 1024710-NLA | 1026192-NLA | 1042989-NLA | 1026761-NLA | 1026768-NLA |
| 19"         | 1024711-NLA | 1026193-NLA | 1042990-NLA | 1026762-NLA | 1026769-NLA |
| 20"         | 1024712-NLA | 1026194-NLA | 1042991     | 1026763-NLA | 1026770-NLA |
| 21"         | 1082654-NLA | 1082952-NLA | 1082953-NLA | 1082954-NLA | 1082955-NLA |
| 22"         | 1082956-NLA | 1082957-NLA | 1082958-NLA | 1082959-NLA | 1082960-NLA |

### #19 - Upright Tube

**Backrests and Hardware**  
**w/o Side Guards or Other Side Guards**  
**#19 - Upright Tube**

|             |              |               |               |               |               |
|-------------|--------------|---------------|---------------|---------------|---------------|
| Back Height | 8 - 11" High | 10 - 14" High | 12 - 16" High | 14 - 18" High | 15 - 19" High |
| Part Number | 1042514-NLA  | 1042513-NLA   | 1043358       | 1042515-NLA   | 1052226-NLA   |

## **Backrests and Hardware**

### **w/o Side Guards or Other Side Guards**

#### **#20 - Backrest Crossmember**

| Chair Width | Part Number | Chair Width | Part Number |
|-------------|-------------|-------------|-------------|
| 12" Wide    | 1042508-NLA | 18" W       | 1050080     |
| 13" Wide    | 1042509-NLA | 19" W       | 1050077     |
| 14" Wide    | 1042510     | 20" W       | 1050074     |
| 15" Wide    | 1042511     | 21" W       | 1085292-NLA |
| 16" Wide    | 1042512     | 22" W       | 1085293-NLA |
| 17" Wide    | 1050083     |             |             |