


Outer Tub Shell and Components

Models 3300, 3600, 3600XL, 3700, 3800, 4300, 5300



SICC_080A

Outer Tub Shell and Components

Models 3300, 3600, 3600XL, 3700, 3800, 4300, 5300

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		PAR7380AG	Back Flow Preventer - Conbraco (1/2") (One Screw in Cap)	1	1
1a		IHPAR7383C	Bonnet	1	1
	D	IHPAR7384C	Kit, Repair	1	1
2		IHPAR7380-AG-W-NLA	NLA - Back-Flow Preventer, "WATTS" 1/2" (Two screws in cap)	1	1
2		IHPAR7382CW-NLA	NLA - Bonnet	1	1
		IHPAR7384CW-NLA	NLA - Kit, Repair	1	1
3		PAR1136	Elbow, 90 Degrees - Chrome	1	1
4		PAR6160D	Washer, Nylon, Backflow Preventer and Chrome Elbow	1	1
5		PAR5080BC	Wand, Shower	1	1
6		PAR5100AC-NLA	NLA - Swivel Connector, Shower Wand	1	1
7		PAR9167B	Holder, Shower Wand	1	1
8		PAR5087B-NLA	NLA - Wand, Disinfectant	1	1
9		PAR9170B	NLA - Holder, Disinfectant Wand	1	1
10	B	PAR9167A	Connector, Chrome Elbow, Hose	1	1
11		IHPAR3645	Leg (6")	1	1
11		IHPAR3646	Leg (10")	1	1
12		PAR9269B	Tub Catch, Chrome	1	1
13		PAR5291-NLA	NLA - Grommet, Rubber	1	1
14		PAR5118G	Hose, Shower and Disinfectant Wand	1	1
15	A	IHPAR5940BC1	NLA - Molding, Trim, Almond	1	1
15	A	IHPAR5940BC	Molding, Trim - White	1	1
15	E	IHPAR1904W	Molding, Trim - White	1	1
16	C	1146438	Preventer, External, Back Flow	1	1
NOTE: A - Used on Base of all freestanding Tubs and end panels on recessed Tubs. B - For Disinfectant Wand. C - For Tubs sold in LA County, Wisconsin, and Oregon only. D - Not Shown in photos. E - Used on 3600, 3600XL, 2300, 3300, 4300, 5300.					