



# Bedend Assembly - Medium Post Headboard

(Models IH297-P-5 and IH297-W-5)

Model IH297-P-5



Model IH297-W-5



## Bedend Assembly - Medium Post Headboard (Models IH297-P-5 and IH297-W-5)

| Item Number                                                                                                                           | Note | Part Number | Description                             | Quantity Per |          |
|---------------------------------------------------------------------------------------------------------------------------------------|------|-------------|-----------------------------------------|--------------|----------|
|                                                                                                                                       |      |             |                                         | Package      | Assembly |
| 1                                                                                                                                     | B    | A024957     | Headboard Assembly                      | 1            | 1        |
| 2                                                                                                                                     | C    | A024979     | Headboard Assembly                      | 1            | 1        |
| 3                                                                                                                                     | D    | IHB001637   | Bumper, Recess, Grey                    | 1            | 2        |
| 4                                                                                                                                     | D    | IHB004149   | Screw, Phillips Pan Head (#8 x 5/8")    | 1            | 2        |
| 5                                                                                                                                     | B    | IHA022678   | Footboard Assembly with Cut Out         | 1            | --       |
| 6                                                                                                                                     | A,B  | A024953     | Footboard Assembly without Cut Out      | 1            | 1        |
|                                                                                                                                       | C    | IHA022721   | Footboard Assembly with Cut Out         | 1            | 1        |
|                                                                                                                                       | C,A  | A024976     | Footboard Assembly without Cut Out      | 1            | 1        |
| 7                                                                                                                                     | D    | IHA016539   | Headboard Assembly, Medium Post - Right | 1            | 1        |
| 8                                                                                                                                     | D    | IHA016538   | Headboard Assembly, Medium Post - Left  | 1            | 1        |
| 9                                                                                                                                     | D    | IHB002194   | T-Nut, Cap (20-1/4")                    | 1            | 4        |
| 10                                                                                                                                    | D    | IHB007478   | Washer, Plated (1/4")                   | 1            | 8        |
| 11                                                                                                                                    | D    | IHB011312   | Bolt, Hex Head (20-1/4 x 3/4")          | 1            | 4        |
|                                                                                                                                       | A,D  | IHB017097   | Cam, Shallow                            | 1            | 6        |
|                                                                                                                                       | A,D  | IHB014947   | Bolt, Hex Head (20-1/4 x 1-1/2")        | 1            | 4        |
| NOTE: A - Not Shown.<br>B - Use with Model IH297-P-5.<br>C - Use with Model IH297-W-5.<br>D - Use with Model IH297-P-5 and IH297-W-5. |      |             |                                         |              |          |