


Back Cane Hardware w/ Lumbar Support

(Fixed Height or Folded Back Canes)



Back Cane Hardware w/ Lumbar Support (Fixed Height or Folded Back Canes)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1172974	Kit, Cane Hardware w/ Lumbar Support (Fixed or Foldable) (Pair)	1	1
2			Bracket, Cane, Rear	1	2
3			Lever, Indexable	1	2
4			Grip, Cane	1	2
5			Insert, Rigidizer	1	2
6			Screw, Cap, Socket Head (M6 x 35mm)	1	2
7			Locknut (M6 x 1)	1	2
8			Screw, Phillips Pan Head (M4.8 X 16mm)	1	2