




# Drawers - Bedside, Chest, Wardrobe, Armoire Torrington



## Drawers - Bedside, Chest, Wardrobe, Armoire Torrington

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1165688	Drawer Assembly, Complete (19-1/8")	1	1
4			Drawer Front Assembly (19-1/8")	1	1
5			Drawer Assembly (19-1/8")	1	1
6			Screw, Phillips Flat Head (#6 x 1")	1	4
2	A	1165689	Drawer Assembly, Complete (28")	1	1
4			Drawer Front Assembly (28")	1	1
5			Drawer Assembly (28")	1	1
6			Screw, Phillips Flat Head (#6 x 1")	1	6
3	A	1165690	Drawer Assembly, Complete (30")	1	1
4			Drawer Front Assembly (30")	1	1
5			Drawer Assembly, Complete (30")	1	1
6			Screw, Phillips Flat Head (#6 x 1")	1	6

NOTE: A - Part includes items 4, 5, and 6.

Furniture Style	Number of Drawers	Part Number
IH28430	Bedside, Single Drawer	1165688
IH28434	Three Drawer Cabinet	1165688
IH28436	Three Drawer Wide Cabinet	1165689
IH28438	Four Drawer Wide Cabinet	1165689
IH28443	Six Drawer Cabinet	1165688
IH28470-2	Two Drawer Wardrobe	1165688
IH28472-2	Two Drawer, Two Door Wardrobe	1165690
IH28490	Armoire, Three Drawers	1165689