




Anti-Tipper Casters and Drive Wheels



Anti-Tipper Casters and Drive Wheels

| Item Number | Note | Part Number | Description | Quantity Per | |
|--|------|-------------|---|--------------|----------|
| | | | | Package | Assembly |
| 1 | A | 1179371-NLA | NLA - Kit, Tire w/ Hardware | 1 | 1 |
| | | | Tire (76 x 318mm) | 1 | 1 |
| | | | Wheel Hardware | 1 | 1 |
| | | | Instruction Sheet, Wheel Installation | 1 | 1 |
| 2 | B | 1179369-NLA | NLA - Kit, Wheel Attaching Hardware | 1 | 1 |
| 3 | | | Hubcap, Drive Wheel | 1 | 1 |
| 4 | | | NLA - Decal, Invacare 3-D Logo (1-19/32" L) | 1 | 1 |
| 5 | | | Locknut (1/2-20) | 1 | 1 |
| 6 | | | Washer (1/2 x 7/8 x 3/64") | 1 | 1 |
| | | | Instruction Sheet, Wheel Installation | 1 | 1 |
| 7 | C | 1179370 | NLA - Kit, Anti-Tip Wheels w/ Hardware | 1 | 1 |
| 8 | | | Wheel, Anti-Tip, Black (3") | 1 | 2 |
| 9 | | | Screw, Socket Button Head Cap (M6 x 1 x 45mm) | 1 | 2 |
| 10 | | | Bushing (1/4 x 3/8 x 1/13") | 1 | 2 |
| 11 | | | Washer (1/4 x 1/2 x 1/16") | 1 | 4 |
| 12 | | | Nut, Lock (M6 x 1mm) | 1 | 2 |
| | | | Instruction Sheet, Anti-Tip Wheels | 1 | 1 |
| NOTE: A - Includes Tire and Wheel Hardware. B - Includes items 3 - 6. C - Includes items 8 - 12. | | | | | |