


# Seat Pan Hardware (20")



## Seat Pan Hardware (20")

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1179375-NLA	NLA - Kit, Seat Pan Hardware (20" W)	1	1
			Outer Shell, Seat Pan - Right	1	1
			Outer Shell, Seat Pan - Left	1	1
			Screw, Button Head Cap (M6 x 40mm)	1	4
			Screw, Round Head Square Neck (M6 x 1 x 16mm)	1	4
			Screw, Pan Head, Tap (8-16 x 1/2")	1	8
			Spacer, Nylon (7 x 18 x 22.4mm)	1	4
			Washer (6.4 x 12.5 x .8)	1	8
			Locknut (M6 x 1)	1	8

NOTE: A - Includes items 2 - 9.