


Anti-Tipper Assembly



Anti-Tipper Assembly

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1159780-NLA	NLA - Kit, Anti-Tip Assembly	1	1
2			NLA - Bracket, Anti-Tip - Right	1	1
3			NLA - Bracket, Anti-Tip - Left	1	1
4			NLA - Anti-Tipper	1	2
5			Screw, Socket Head (1/4-28 x 7/8")	1	4
6			Spacer, Crossbrace (1/4 x 7/16 x 5/16)	1	4

NOTE: A - Includes items 2, 3, 4, and 5.

Anti-Tip Rear Seat To Floor Reference Chart

Rear Seat to Floor	CRF22SP, CRF6722 (22" Wheel)	ZEROG24, CRF6724 SPLX24, FN1624 (24" Wheel)	ZEROG25, CRF6725 SPLX25, FN1625 ZERO26, CRF6726, FN1626 (25 & 26" Wheels)
14"	1159780	1159780	1159780
14-1/2"	1159780	1159780	1159780
15"	1159780	1159780	1159780
15-1/2"	1159780	1159780	1159780
16"	1159780	1159780	1159780
16-1/2"	1159781-NLA	1159780	1159780
17"	1159781-NLA	1159780	1159780
17-1/2"	1159781-NLA	1159781-NLA	1159780
18"	1159781-NLA	1159781-NLA	1159780
18-1/2"	1159781-NLA	1159781-NLA	1159781-NLA
19"	1159781-NLA	1159781-NLA	1159781-NLA
19-1/2"	1159781-NLA	1159781-NLA	1159781-NLA
20"	1159781-NLA	1159781-NLA	1159781-NLA