

Adjustable Height w/o Push Handles (ABH)

For Back Brace and Indexable Lever, see category titled "Upholstery - Back Brace and Hardware" page.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1172977	Kit, Cane Hardware w/o Push Handle (Adjustable Height)	1	1
2			Bracket, Cane, Rear	1	2
3			Bushing, Flanged	1	2
4			Button, Plug (22mm)	1	2
5			Screw, Socket Head Cap (M6 x 1 x 30)	1	2
6			Screw, Cap, Socket Head (M6 x 35mm)	1	2
7			NLA - Locknut (M6 x 1)	1	2
8	B	1172976	Kit, Cane Bracket Attaching Hardware (Pair)	1	1
7			NLA - Locknut (M6 x 1)	1	4
9			Screw, Socket Button Head (M6 x 53mm)	1	2
10			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
11		1169322	Insert, Back Cane, Adjustable Height	1	2
12		1172983	Back Cane, Adjustable Height, Straight (12 - 16")	1	2
12		1172985	Back Cane, Adjustable Height, Straight (15 - 19")	1	2
12		1172996	Back Cane, Adjustable Height, 8 Degree (12 - 16")	1	2
12		1172998	Back Cane, Adjustable Height, 8 Degree (15 - 19")	1	2

NOTE: A - Kit includes items 2 - 7.
B - Kit includes items 7, 9, and 10.

Back Canes - Adjustable Height w/o Push Handles (ABH)

Back Cane Style	Item #1 Hardware Part No.	Item #8 Cane Attaching Hardware Part No.	Item #12 Back Cane Part No.
ABH1216 (8 Degree)	1172977	1172976	1172996
ABH1216 (Straight)	1172977	1172976	1172983
ABH1519 (8 Degree)	1172977	1172976	1172998
ABH1519 (Straight)	1172977	1172976	1172985
ABH1620 (8 Degree)	1172977	1172976	N/A
ABH1620 (Straight)	1172977	1172976	N/A