


Amigo Rollator



Amigo Rollator

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	13571	Kit, Wheel with Fork Assembly SL200-Narrow - Front	1	1
		80079	Wheel, Front	1	2
		72210	Sleeve	1	2
		62007	Washer	1	2
		61016	Locknut (M8)	1	2
		60040	Screw	1	2
		13550	Fork	1	2
		70050-NLA	NLA - Cap, Black Plastic - Front Fork	1	2
2		13573	Wheel Assembly with Hardware, SL200-Narrow - Front	1	2
		12521	Wheel, SL200-Narrow - Front	1	2
3		70050-NLA	NLA - Cap, Black Plastic - Front Fork	1	2
4	B	13574	Kit, Brake Sleeve	1	1
		70308	Sleeve, Brake	1	2
		62066	Block, Brake	1	2
		60099	Screw, Oval Head (M2.5 x .45 x 10mm)	1	2
5	C	12533	Kit, Wide Wheel Assembly with Hardware - Rear	1	1
		80071	Wheel, Wide - Rear	1	2
		60041	Screw (M8 x 80mm)	1	2
		60148	Screw, Stop (M8 x 80mm)	1	2
		61022	Nut	1	2
		62007	Washer	1	2
		63040	Sleeve	1	2
		70066	Bushing	1	2
		70093	Bushing	1	2
6		12660	Shoe, Brake	1	2
7		13040	Kit, Handle DABS - Left	1	1
7		13039	Kit, Handle DABS - Right	1	1
8		13464-L	Brake Handle Assembly - Left	1	1
8		13464-R	Brake Handle Assembly - Right	1	1
9		80078	Wire (1115mm)	1	2
10		80077	Wire Cover (1050mm)	1	2
11		12509	Grip Anatomic with Cover (Pair)	2	1
12		12677	Grip, Soft (220mm) (Pair)	2	1
13		12545	Seat Assembly with Hardware, Flip-up	1	1
14		12656	Lock, Folding	1	1
		70202	Latch, Safety w/o Holder	1	1
15		12615	Latch, Safety, Seat	1	1
16		12622	Kit, Folding Handle	1	1
17		12665	Adjustment Knob and Screw (Pair)	2	1
18	D	12648	Folding Hinge Assembly (Pair)	2	1
		70179	Spacer, Tube	1	1
		13337-16	Stopper, Seat	2	2
		60001	Screw (M8 x 80mm)	2	2
		60109	Screw, Seat Support (M5 x .8 x 30mm)	2	2
19		12647	Rail, Folding with Hardware	1	1

Amigo Rollator

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
		82057	Instruction Sheet	1	1
NOTE: A - Includes: Front Wheel, Sleeve, Washer, Fork, Screw, Nut and Cap. B - Includes: Brake Sleeve, Brake Block and Screw. C - Includes: Rear Wheel, Screw, Stop Screw, Nut, Washer, Sleeve and Bushings. D - Includes: Seat Stopper, Spacer and Screws.					