


Seat & Arms

Leo 3W, Leo 4W



Seat & Arms Leo 3W, Leo 4W

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1173792-NLA	NLA - Arm Assembly - Right	1	1
1	A	1173793-NLA	NLA - Arm Assembly - Left	1	1
1	B	1506314	Arm Assembly - Right	1	1
1	B	1506313	Arm Assembly - Left	1	1
2		1164833	Seat Slide Assembly	1	1
3		1164835-NLA	NLA - Seat and Back Assembly	1	1
3	E	1578946	Seat and Back Assembly w/ Double White Stitching	--	1
4		1164834-NLA	NLA - Seat Bracket Assembly	--	1
5		1165106	Post, Seat Guide w/ Hardware	--	1
6		1173789	Armrest	--	2
7	C,F	1168463-NLA	NLA - Cover, Back, Leo - Black Only	--	1
8	C,G	1168464-NLA	NLA - Cover, Seat, Leo - Black Only	--	1
	D	1507005	Knob, Tightening	--	1

NOTE: A - For Three or Four Wheel Scooters manufactured before 08/08/2011.
 B - For Three or Four Wheel Scooters manufactured after 08/07/2011.
 C - For complete Seat Assembly, see item 3 (1164835).
 D - Not shown in the parts illustration.
 E - After February 2014, Serial Number 14CT000295.
 F - Before Serial Number 14CT000294.
 G - Before Serial Number 14CT000294.