


# Motor & Transaxle



## Motor & Transaxle

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1164829	Core, EMI (Pair)	2	1
2	A	1506318	Lever, Brake w/ Hardware	1	1
3	B	1164828	Motor w/ Brake and Hardware	1	1
4	C	1164827	Transaxle w/ Hardware	1	1
	D	1534599	Carbon Brush	1	2
NOTE: A - Includes (2) screws and washers. B - Includes bolt and guide. C - Includes (2) key shafts. D - Item is not shown in the illustration.					