



Back Upholstery and Hardware

9000SL / 9000XT / 9000XDT


If a wheelchair was originally ordered with the RB, (reinforced back), please refer to quote SM-100924021021 to order a replacement back upholstery.




9000SL
Back Upholstery
Wrap-Around
15-19" - Nylon/Vinyl




9000SL
Back Upholstery
Screw-On
15-19" - Nylon




9000SL
Back Upholstery
Screw-On
15-19" - Vinyl




9000XT
Back Upholstery
Wrap-Around
15-19" - Nylon/Vinyl




9000XT
Back Upholstery
Screw-On
15-19" - Vinyl



9000XT
Back Upholstery
Screw-On
15-19" - Vinyl



9000XDT
Back Upholstery
Wrap-Around
15-19" - Nylon/Vinyl



9000XDT
Back Upholstery
Screw-On
15-19" - Vinyl

Back Upholstery and Hardware

9000SL / 9000XT / 9000XDT

If a wheelchair was originally ordered with the RB, (reinforced back), please refer to quote SM-100924021021 to order a replacement back upholstery.

9000SL Nylon Screw-On Back Upholstery

Nylon - Before 7/12/99

Chair Width	15" Back Height	16" Back Height	17" Back Height	18" Back Height	19" Back Height
16" W	8881057251 See Note(s) A	8881057252-NLA See Note(s) A	8881057253-NLA See Note(s) A	8881057254-NLA See Note(s) A	8881057255-NLA See Note(s) A
18" W	8881057256-NLA See Note(s) A	8881057257-NLA See Note(s) A	8881057258-NLA See Note(s) A	8881057259-NLA See Note(s) A	8881057260-NLA See Note(s) A
20" W	N/A	N/A	N/A	8881066271-NLA See Note(s) A	8881066272-NLA See Note(s) A

NOTE: A - Only U240 Black Nylon Upholstery is available.

NOTE: Hardware is not included. See chart titled, "Hardware for Back Upholstery - 9000SL."

9000SL Nylon Screw-On Back Upholstery

Nylon - After 7/11/99

Chair Width	Back Height 15 - 17"	Back Height 17 - 19"
16" W	1127480	1127481
18" W	1127482	1127483
20" W	1127484	1127485

NOTE: Upholstery color is black.

NOTE: Hardware is included

9000SL Screw-On Back Upholstery

Vinyl - After 7/11/99

Chair Width	Back Height 15 - 17"	Back Height 17 - 19"
16" W	8881127500	8881127501
18" W	8881127502	8881127503
20" W	8881127504	8881127505

NOTE: When ordering, specify upholstery color, (U67).

NOTE: Hardware is included.

Back Upholstery and Hardware

9000SL / 9000XT / 9000XDT

If a wheelchair was originally ordered with the RB, (reinforced back), please refer to quote SM-100924021021 to order a replacement back upholstery.

9000SL Wrap-Around Upholstery

(For ANGLE BACK TYPE Option)

Chair Width	Back Height 15-19" - Nylon	Back Height 15-19" - Vinyl
16" W	1130878 See Note(s) A	8881130886 See Note(s) B
18" W	8881130880 See Note(s) A	8881130888 See Note(s) B
20" W	1130882 See Note(s) A	8881130890 See Note(s) B

NOTE: A - Upholstery color is Black.
B - Specify Upholstery Color (U67), when ordering.

NOTE: Hardware is included.

Hardware for Back Upholstery - 9000SL

Part Number	Description	Quantity
1041541 See Note(s) A	Screws w/ Washers	1 package of 12
1018239 See Note(s) B	Tie Wrap	2

NOTE: A - Hardware is not sold individually.
B - Sold individually; assembly requires (2).

Back Upholstery and Hardware

9000SL / 9000XT / 9000XDT

If a wheelchair was originally ordered with the RB, (reinforced back), please refer to quote SM-100924021021 to order a replacement back upholstery.

9000XT & 9000XDT Back Upholstery

Chair Width	Nylon Wrap-Around Upholstery 15 - 19" Back Height	Vinyl Wrap-Around Upholstery 15 - 19" Back Height	Vinyl Screw-On Upholstery 15 - 19" Back Height
14" W	1130876 See Note(s) A	8881130884 See Note(s) B	8881130892 See Note(s) B
15" W	1130877 See Note(s) A	8881130885 See Note(s) B	8881130893 See Note(s) B
16" W	1130878 See Note(s) A	8881130886 See Note(s) B	8881130894 See Note(s) B
17" W	1130879 See Note(s) A	8881130887 See Note(s) B	8881130895 See Note(s) B
18" W	8881130880 See Note(s) A	8881130888 See Note(s) B	8881130896 See Note(s) B
19" W	1130881 See Note(s) A	8881130889 See Note(s) B	8881130897 See Note(s) B
20" W	1130882 See Note(s) A	8881130890 See Note(s) B	8881130898 See Note(s) B
22" W	1130883 See Note(s) A	8881130891 See Note(s) B	8881130899 See Note(s) B

NOTE: A - Upholstery color is black.
B - When ordreing, specify upholstery color, (U67).

NOTE: Hardware is included.

Hardware for Back Upholstery - 9000XT

Part Number	Description	Quantity
1041541 See Note(s) A	Screws w/ Washers	1 package of 12
1018239 See Note(s) B	Tie Wrap	2

NOTE: A - Hardware is not sold individually.
B - Sold individually; assembly requires 2.

Hardware for Back Upholstery - 9000XDT

Part Number	Description	Quantity
1041541 See Note(s) A	Screws w/washers	1 package of 12

NOTE: A - Hardware is not sold individually.