

Manual Lift - Base

NOTE: All fasteners must be torqued to IFI/DIN standards.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1156139	Leg Assembly	1	2
2		1175245-NLA	NLA - Cover, Leg Protection	1	2
3		1156233-NLA	NLA - Shroud, Manual Base	1	1
4		1023632	Clip, Shroud	1	3
5	A	1177622-NLA	NLA - Kit, Rear Caster w/ Hardware, Locking (Pair)	1	1
12			Screw, Socket, Countersunk Head (M10 x 1.5 x 85mm)	1	2
13			Spacer (41.3 x 10.2 x 15mm)	1	2
14			Caster, Rear Locking	1	2
23			Locknut (M10 x 1.5mm)	1	2
6	B	1177623-NLA	NLA - Kit, Front Caster w /Hardware (Pair)	1	1
15			Screw, Socket, Countersunk Head (M10 x 1.5 x 70mm)	1	2
16			Spacer (40 x 10.5 x 3mm)	1	2
17			Caster, Front	1	2
23			Locknut (M10 x 1.5mm)	1	2
7	C	1177624-NLA	NLA - Kit, Leg to Crossmember Mounting Hardware (Pair)	1	1
18			Bushing, Bronze	1	4
19			Bushing, Base and Leg	1	2
20			Screw, Socket Button Head, (M10 x 1.5 x 22mm)	1	2
31			Lock Washer, Split (M10)	1	2
8	D	1177626-NLA	NLA - Kit, Index Plate w /Hardware	1	1
21			Washer (10.5 x 20 x 2mm)	1	2
22			Index Plate Assembly	1	1
23			Locknut (M10 x 1.5mm)	1	2
9	E	1177627-NLA	NLA - Kit, Handle Mount	1	1
24			Handle, Mount	1	1
25			Spring, Compression (1 x 1-1/16 x 1/8")	1	1
26			Screw, Shoulder Hex Socket Head w/ Patch (1/2" x 1-1/4", 3/8-16)	1	1
10	F	1177629-NLA	NLA - Kit, Linkage (Pair)	1	1
27			Screw, Socket Button Head w/ Patch (M8 x 1.25 x 25mm)	1	4
28			Tie Rod Assembly	1	2
11	G	1177628	NLA - Kit, Handle w /Grip	1	1
29			Handle	1	1
30			Grip, Cane	1	1

NOTE: A - Includes items 12-14 and 23.
 B - Includes items 15-17 and 23.
 C - Includes items 18-21.
 D - Includes items 21-23.
 E - Includes items 24-26.
 F - Includes items 27-28.
 G - Includes items 29-30.