


Back Canes - Adjustable Angle w/ Hardware (20" H)



Back Canes - Adjustable Angle w/ Hardware (20" H)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1169696	Kit, Back Cane Hardware, Adjustable Angle (20") (Pair)	1	1
2			Plate, Adjustable Angle, Outer - Right	1	1
2			Plate, Adjustable Angle, Outer - Left	1	1
3			Plate, Adjustable Angle, Inner - Right	1	1
3			Plate, Adjustable Angle, Inner - Left	1	1
4			Backrest, Excenter	1	2
5			Support, Seat Tube, Rear	1	2
6			Spacer, Coved (1/4 x 1 x 13/32")	1	8
7			Insert, Back Cane, Adjustable Height (20" H)	1	2
8			NLA - Screw, Socket Head Cap (M6 x 1 x 40mm)	1	4
9			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
10			Screw, Socket Button Head (M6 x 53mm)	1	2
11			NLA - Locknut (M6 x 1)	1	2
12			Washer (1/4 x 9/16 x 1/16")	1	2
13			Screw, Socket Head Cap (M6 x 1 x 30)	1	2
3		1179667	Plate, Adjustable Angle, Inner - Right	1	1
3		1179668	Plate, Adjustable Angle, Inner - Left	1	1
5		1533913	Support, Seat Tube, Rear	1	2
7		1169635	Insert, Back Cane, Adjustable Height (20" H)	1	2
NOTE: A - Includes items 2 - 13.					