



Shower Chair - Compact Blow Molded Fixed Height

Model 50-1, 3
Without Back



Model 60-1, 3
With Back



Shower Chair - Compact Blow Molded Fixed Height

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1079341-NLA	NLA - Seat	1	1
2	C	1079340-NLA	NLA - Back	1	1
3		1079313-NLA	NLA - Leg, Top	1	1
4		1079312-NLA	NLA - Leg, Bottom	1	1
5	C	1083245-NLA	NLA - Tube, Backrest Support	1	2
6		40186X030	Tip, Rubber, Grey (1")	1	4
7		1027096	Plug Button, Grey (7/8")	1	2
	A,B,D	1074249-NLA	NLA - Kit, Hardware	1	1
8			NLA - Nut, Wing (1/4-20)	1	1
9			NLA - Washer, External Tooth Lock (1-1/4")	1	1
10			Screw, Phillips Pan Head (1/4-20 x 1-3/8")	1	1
11			NLA - Screw, Slotted Hex Washer Head (1/4-14 x 1-1/2")	1	4
12			NLA - Washer, Steel Flat (1/4 x 5/8 x 1/16")	1	1
	B,C,E	1079912-NLA	NLA - Kit, Hardware	1	1
13			NLA - Nut, Wing (1/4-20)	1	1
14			NLA - Washer, External Tooth Lock (1-1/4")	1	1
15			Screw, Phillips Pan Head (1/4-20 x 1-3/8")	1	1
16			NLA - Screw, Phillips Pan Head (1/4-14 x 2-1/4")	1	4
17			NLA - Washer, Steel Flat (1/4 x 5/8 x 1/16")	1	1
18			NLA - Screw, Phillips Pan Head Tap (#8-18 x 1/2")	1	2
NOTE: A - To be used with Model 50-1,3. B - Not Shown. C - To be used with Model 60-1,3. D - Kit Includes Items 8-12. E - Kit Includes Items 13-18.					