


Footrest and Hardware



Footrest and Hardware

For Footrest Accessories, Refer To Miscellaneous Accessories Category

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1170247-NLA	NLA - Kit, Hardware, Footplate	1	1
3			Clamp, Footrest	1	1
5			Screw, Socket, Flat Head (1/4-28 x 1")	1	2
7			Locknut (1/4-28)	1	2
2	B	1170248-NLA	NLA - Kit, Hardware, Footplate, 2" Spacer Riser	1	1
3			Screw, Socket, Flat Head (1/4-28 x 1")	1	2
4			Spacer (16/64 x 3/4 x 2")	1	2
5			NLA - Screw, Socket Flat Head (1/4-28 x 1")	1	2
6			Clamp, Footrest	1	1
8	C	1143290-NLA	NLA - Footplate (6-1/2" W)	1	1
8	D	1143291-NLA	NLA - Footplate (7-1/2" W)	1	1
8	E	1143292-NLA	NLA - Footplate (8-1/2" W)	1	1
8	F	1143293	Footplate (9-1/2" W)	1	1
8	G	1143286-NLA	NLA - Footplate (10-1/2" W)	1	1
9	J	8881143283	Tube, Support, Footrest, Tapered, Black	1	2
9		1143282-NLA	NLA - Tube, Support, Footrest, V Front, Black	1	2
10	H,J	8881146354-NL A	NLA - Footrest, Tubular w/ Holes, Black	1	1
10	H,J I	8881143278 See Chart	Footrest, Tubular w/o Holes, Black Tubular Footrest Covers	1 1	1 1
NOTE: A - Includes items 3, 5, and 7. B - Includes items 3, 4, 5, and 6. C - For use on chairs 12 - 13" wide. D - For use on chairs 14 - 15" wide. E - For use on chairs 16 - 17" wide. F - For use on chairs 18 - 19" wide. G - For use on chairs 20" wide. H - Use with Tubular Footrest Options, See below chart. I - Not shown in illustration. J - Order as a "special" order. Reference quote #090723085034 when ordering. The Serial Number and chair width must be added to order for the replacement service part also.					

Tubular Footrest Covers

	12 - 14" W	15" W	16" W	17" W	18" W	19" W	20" W
Plastic Footrest Cover	1145654-NLA	1145655-NLA	1145656-NLA	1145657-NLA	1145658-NLA	1145659-NLA	1145660-NLA
Carbon Fiber Footrest Cover	1145647-NLA	1145648-NLA	1145649-NLA	1145650-NLA	1145651-NLA	1145652-NLA	1145653-NLA