


Rollator (Four Wheeled) - Model 65650 and 65650R



The photo represents a previous model of the 65650 Rollator, Review all Footnotes for additional information before ordering any part.

Rollator (Four Wheeled) - Model 65650 and 65650R

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1142628	Grip, Hand (Pair)	1	1
2	B,F	1142629-NLA	NLA - Brake Handle w/ Hardware Assembly	1	2
3	C,G	1142630-NLA	NLA - Brake Assembly	1	2
4		1142632-NLA	NLA - Back Support Assembly, Curved	1	1
5		1142633-NLA	NLA - Seat w/ Hardware	1	1
6		1142638-NLA	NLA - Basket	1	1
7		1142639-NLA	NLA - Connection Linkage Assembly	1	2
8		1142637-NLA	NLA - Package, Strap	10	1
9	D	1142636-NLA	NLA - Kit, Wheel w/ Hardware - Rear	1	2
10	F	1142634-NLA	NLA - Fork Assembly, Front	1	2
11		1142635-NLA	NLA - Kit, Wheel w/ Hardware - Front	1	2
12	E	1144495-NLA	NLA - Kit, Adjustable Knob (Pair)	2	1
12		1169007	Kit, Adjustable Knob (Pair)	2	1

NOTE: A - Includes (2) Grips.
 B - Includes Handle with attaching Hardware.
 C - Includes Handle with Cable, Fork, Brake Pad and Hardware.
 D - Includes Rear Wheel with Hardware.
 E - Part is no longer available. See #1169007 for replacement.
 F - Replacement wheel has different amount of Spokes that does not impact performance. It is recommended to replace Fork Assembly in pairs to maintain appearance.
 G - New Handgrips must be ordered as existing Handgrips must be cut off to install new Brake Assembly.