


Rails



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1557047	Rail, Aluminum, Midi (3 m / 9 ft., 10 in.)	--	--
1		1557052	Rail, Aluminum, Midi (3.5 m / 11 ft., 6 in.)	--	--
1		1557072	Rail, Aluminum, Midi (4 m / 13 ft., 2 in.)	--	--

Rails Mounting Kits

Part Number	Qty.	Capacity	Mounting Used With
1557047	1	205 kg. (450 lbs.)	3M (9 ft., 10 in.) Rail, Wall, Ceiling, and Upright
1557047	3	205 kg. (450 lbs.)	3M (9 ft., 10 in.) Rail H-Frame Wall, Ceiling, and Upright
1557047	1	250 kg. (550 lbs.)	3M (9 ft., 10 in.) Rail Wall, Ceiling and Upright
1557047	3	250 kg. (550 lbs.)	3M (9 ft., 10 in.) Rail H-Frame Wall, Ceiling, and Upright
1557052	1	250 kg. (550 lbs.)	3.5M (11 ft., in.) Rail, Wall and Upright

Rails

Rails Mounting Kits

Part Number	Qty.	Capacity	Mounting Used With
1557052	2	250 kg. (550 lbs.)	3.5 (11 ft., 6 in.) x 4M (13 ft., 2 in.) Room Covering Wall and Upright
1557072	1	205 kg. (450 lbs.)	4M (13 ft., 2 in.) Rail Wall, Ceiling, and Upright
1557072	3	205 kg. (450 lbs.)	4M (13 ft., 2 in.) Room Covering, Wall, Ceiling and Upright
1557072	1	250 kg. (550 lbs.)	4M (13 ft., 2 in.) Rail, Ceiling
1557072	1	250 kg. (550 lbs.)	3.5 (11 ft., 6 in.) x 4M (13 ft., 2 in.) Room Covering Wall and Upright
1557072	3	250 kg. (550 lbs.)	4M (13 ft., 2 in.) Room Covering, Ceiling