


Mounting Brackets



| Item Number | Note | Part Number | Description | Quantity Per | |
|-------------|------|-------------|---|--------------|----------|
| | | | | Package | Assembly |
| 1 | | 1557046 | Package, End Stops and Wall Mount Brackets | 1 | 1 |
| 2 | | 1557109 | Package, End Stops and Ceiling Mount Brackets | 1 | 1 |
| 3 | | 1557111 | Package, End Stops and Room Covering Wall Brackets | 1 | 1 |
| 4 | | 1557112 | Package, End Stops and Room Covering Ceiling Brackets | 1 | 1 |
| 5 | | 1557110 | Package, End Stops and Upright Support Mounts | 1 | 1 |
| 6 | | 1557113 | Package, End Stops and Room Covering Upright Mounts | 1 | 1 |

Mounting Brackets

Brackets Used with Mountings

| Part Number | Capacity | Mounting Description | Mounting Part No. |
|-------------|--------------------|--|-------------------|
| 1557046 | 250 kg. (550 lbs.) | 3.5M (11 ft., 6 in.) Ceiling | 1557022 |
| 1557046 | 250 kg. (550 lbs.) | 3M (9 ft., 10 in.) Wall | 1556975 |
| 1557046 | 205 kg. (450 lbs.) | 4M (13 ft., 2 in.) Wall | 1556923 |
| 1557046 | 205 kg. (450 lbs.) | 3M (9 ft., 10 in.) Wall | 1556861 |
| 1557109 | 250 kg. (550 lbs.) | 4M (13 ft., 2 in.) Ceiling | 1557023 |
| 1557109 | 250 kg. (550 lbs.) | 3M (9 ft., 10 in.) Ceiling | 1556987 |
| 1557109 | 205 kg. (450 lbs.) | 4M (13 ft., 2 in.) Ceiling | 1556924 |
| 1557109 | 205 kg. (450 lbs.) | 3M (9 ft., 10 in.) Ceiling | 1556880 |
| 1557110 | 250 kg. (550 lbs.) | 3.5M (11 ft., 6 in.) Upright | 1557042 |
| 1557110 | 250 kg. (550 lbs.) | 3M (9 ft., 10 in.) Upright | 1557017 |
| 1557110 | 205 kg. (450 lbs.) | 4M (13 ft., 2 in.) Upright | 1556962 |
| 1557110 | 205 kg. (450 lbs.) | 3M (9 ft., 10 in.) Upright | 1556900 |
| 1557111 | 250 kg. (550 lbs.) | 3.5 (11 ft., 6 in.) x 4M (13 ft., 2 in.) Room Covering, Wall | 1557043 |
| 1557111 | 250 kg. (550 lbs.) | 3M (9 ft., 10 in.) Room Covering, Wall | 1557018 |
| 1557111 | 205 kg. (450 lbs.) | 4M (13 ft., 2 in.) Room Covering, Wall | 1556963 |
| 1557111 | 205 kg. (450 lbs.) | 3M (9 ft., 10 in.) Room Covering, Wall | 1556901 |
| 1557112 | 250 kg. (550 lbs.) | 4M (13 ft., 2 in.) Room Covering, Ceiling | 1557044 |
| 1557112 | 250 kg. (550 lbs.) | 3M (9 ft., 10 in.) Room Covering, Ceiling | 1557019 |
| 1557112 | 205 kg. (450 lbs.) | 4M (13 ft., 2 in.) Room Covering, Ceiling | 1556973 |
| 1557112 | 205 kg. (450 lbs.) | 3M (9 ft., 10 in.) Room Covering, Ceiling | 1556902 |
| 1557113 | 250 kg. (550 lbs.) | 3M (9 ft., 10 in.) Room Covering, Upright | 1557020 |
| 1557113 | 250 kg. (550 lbs.) | 3.5 (11 ft., 6 in.) x 4M (13 ft., 2 in.) Room Covering, Upright | 1556845 |
| 1557113 | 205 kg. (450 lbs.) | 3M (9 ft., 10 in.) Room Covering, Upright | 1556903 |
| 1557113 | 205 kg. (450 lbs.) | 4M (13 ft., 2 in.) Room Covering Upright | 1556974 |