




# Hill Holder Wheel Locks

## 9695, Conventional Arm, 9XT/9SL



## Hill Holder Wheel Locks

### 9695, Conventional Arm, 9XT/9SL

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1034210	Wheel Lock Assembly, Hill Holder, Push to Lock - Right	1	1
1	A	1034211	Wheel Lock Assembly, Hill Holder, Push to Lock - Left	1	1
1	B	1025429-NLA	NLA - Wheel Lock Assembly, Hill Holder, Pull to Lock - Right	1	1
1	B	1025430-NLA	NLA - Wheel Lock Assembly, Hill Holder, Pull to Lock - Left	1	1
2		47376X024	Tip, Rubber Wheel Lock Handle, Black	1	1
3		1021402	Bolt, Carriage (1/4-20 x 1-1/4")	1	1
4		1021403	Nut, Hex Flang (1/4-20)	1	1
NOTE: A - For use with Adult style chairs B - For use with Hemi style chairs.					
NOTE: Hill holders right and left can only be used on 22" and 24" Pneumatic Tires.					