


Shock Absorbers



Shock Absorbers

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,C	1041419-NLA	NLA - Shock Assembly, (RWD) Standard	1	1
1	A,D	1039660	Shock Assembly, (RWD) Heavy Duty	1	1
1	A,B,E	1070371-NLA	NLA - Shock Assembly, (MWD) Heavy Duty	1	1
2		1041431	Shock Eyelet Assembly	1	1
3		1054481-NLA	NLA - Retainer, Shock Spring, Upper	1	1
4	C	1043595-NLA	NLA - Spring, Shock, Standard	1	1
4	D	1043596-NLA	NLA - Spring, Shock, Heavy Duty	1	1
5		1042466-NLA	NLA - Retainer, Shock Spring, Lower	1	1
6		1039659-NLA	NLA - Absorber, Shock	1	1

NOTE: A - Includes items 2 - 6.
 B - Due to spring compression, the heavy duty shock absorber is not saleable as individual components.
 C - Junior, 200 pound capacity
 D - 250 pound capacity
 E - 400 pound capacity