



Shower Chair - Large Blow Molded Adjustable Height

Model 9680, 9681P
w/o Back
(Discontinued 12/1/04)



Model 9680P, 9681
w/ Back
(Discontinued 12/1/04)



Shower Chair - Large Blow Molded Adjustable Height

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,C	1060375-NLA	NLA - Back	1	1
2		1060376-NLA	NLA - Seat	1	1
3	F	1061419-NLA	NLA - Leg, Bottom	1	1
3	G	1079315-NLA	NLA - Leg, Bottom	1	1
4	F	1061418-NLA	NLA - Leg, Top	1	1
4	G	1079314-NLA	NLA - Leg, Top	1	1
5	E	1065346-NLA	NLA - Leg Extension Assembly	1	4
6		1063124	Anti-Rattle, Grey (1-1/8")	1	4
7		1057797	Tip, Rubber, Light Grey (1-1/8")	1	4
8		1094945	Button, Snap	1	4
9	A,D,J	1044142	Screw, Phillips Pan Head (1/4-20 x 1-3/8")	1	1
10		1008882	Locknut (1/4-20)	1	1
13		1021944-NLA	NLA - Screw, Slotted Hex Washer Head (1/4-14 x 1-1/2")	1	4
		1074248-NLA	NLA - Kit, Hardware	1	1
9	B,C,K	1074249-NLA	Screw, Phillips Pan Head (1/4-20 x 1-3/8")	1	1
12			Locknut (1/4-20)	1	1
13			NLA - Screw, Slotted Hex Washer Head (1/4-14 x 1-1/2")	1	4
			NLA - Kit, Hardware	1	1
9	H	1068492-NLA	Screw, Phillips Pan Head (1/4-20 x 1-3/8")	1	1
10			NLA - Washer, Steel Flat (1/4 x 5/8 x 1/16")	1	1
11			NLA - Washer, External Tooth Lock (1-1/4")	1	1
12			NLA - Nut, Wing (1/4-20)	1	1
13	H,I	1064507-NLA	NLA - Screw, Slotted Hex Washer Head (1/4-14 x 1-1/2")	1	4
			NLA - Back With Hardware	1	1
			NLA - Plug, End	1	2

NOTE: A - To be used with Model 9681.
 B - To be used with Model 9681P.
 C - To be used with Model 9680P.
 D - To be used with Model 9680.
 E - Includes Items 6 and 7.
 F - Before 11/2/98.
 G - After 11/1/98.
 H - Not Shown.
 I - Used with Back before 5/15/97.
 J - Kit Includes Items 9, 12 and 13.
 K - Kit Includes Items 9, 10, 11, 12 and 13.