

Warranty Replacement Parts

Part Number	Description	Quantity per Assembly
1108720E-NLA See Note(s) A	WARR RPLMT - LXN PCB Assembly	1
1108731E-NLA See Note(s) A	WARR RPLMT - ISBN PCB Assembly	1
1106690E-NLA See Note(s) B	WARR RPLMT - Compressor Assembly, Thomas (5LX Plat.)	1
1128741E-NLA See Note(s) C	WARR RPLMT - Compressor (10LX/10LXO2)	1
1130666E See Note(s) G, C	WARR RPLMT - Compressor Assembly (10LX/10LXO2)	1
1143486E See Note(s) D	WARR RPLMT - PCB Assembly (PO2/PO2W)	1
1128742E See Note(s) E	WARR RPLMT - PCB Control Assembly (10LXO2)	1
1075071E-NLA See Note(s) F	WARR RPLMT - Compressor Assembly (6LX/6LXO2)	1
<p>NOTE: A - See "Control Panel and PC Boards (Platinum 5 Models - After 12/5/02)" Category for more information. B - See "Compressor and Mtg. Hardware (Platinum 5 Models)" Category for more information. C - See "Compressor and Mtg. Hardware (10 Liter Models)" Category for more information. D - See "Control Panel and PC Boards (Perfecto2/PerfectoW)" Category for more information. E - See "Control Panel - Rear View" Category for more information. F - See "Compressor and Mtg. Hardware (Non Platinum 5 Models)" Category for more information. G - When ordering 1130666E, you must order Kit part number 1186230.</p> <p>When ordering 1130666E, you must order Kit part number 1186230.</p>		