


Versair Adjustable Tension Back Upholstery Assemblies and Privacy Flaps



Front view of Back Assembly



Back Adjuster shown for Reference only

Privacy Flap shown attached to Adjuster



Wheelchair Shown for Reference only

Versair Adjustable Tension Back

Versair Adjustable Tension Back Upholstery Assemblies and Privacy Flaps

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1169271	Upholstery Assembly, w/ Hardware, Back, Versair, Adj. Tension Black Mesh (16-21" W x 15-17" H)	1	1
1	A	1169272	Upholstery Assembly, w/ Hardware, Back, Versair, Adj. Tension Black Mesh (21-26" W x 15-17" H)	1	1
1	A	1169273	Upholstery Assembly, w/ Hardware, Back, Versair, Adj. Tension Black Mesh (16-21" W x 17-19" H)	1	1
1	A	1169274	Upholstery Assembly, w/ Hardware, Back, Versair, Adj. Tension Black Mesh (21-26" W x 17-19" H)	1	1
2		1179453	Upholstery, Privacy Flap, Versair, Black Nylon, Standard (16 - 21" W)	1	1
2		1179454	Upholstery, Privacy Flap, Versair, Black Nylon, Large (21 - 26" W)	1	1

NOTE: A - Assembly includes Back Upholstery, Right and Left Adjusters, Hardware and Instruction Sheet. Available as a Back Upholstery Assembly only. Individual components are not serviceable.