

Veranda Wheelchair

Footrest packaged with the wheelchair match frame finish. See Section Titled, "Veranda Only" for T93H/T94H Series Front Rigging replacement Options.

At this time, Upholstery replacements for 16" and 20" Chairs are unavailable.

The new front caster used on a 16 and 20" Chair was changed as of March 8, 2011 and for the 18" Chair was changed as of June 21, 2011. "New style caster is available through The Aftermarket Group as TAGRP265093"



Veranda Wheelchair

Footrest packaged with the wheelchair match frame finish. See Section Titled, "Veranda Only" for T93H/T94H Series Front Rigging replacement Options.

At this time, Upholstery replacements for 16" and 20" Chairs are unavailable.

The new front caster used on a 16 and 20" Chair was changed as of March 8, 2011 and for the 18" Chair was changed as of June 21, 2011. "New style caster is available through The Aftermarket Group as TAGRP265093"

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1132515	Kit, Hand Grip (Pair)	1	1
2		1152605	Kit, Back Upholstery w/ Hardware, Removable Arm (18 x 16")	1	1
2		1157814	Kit, Back Upholstery w/ Hardware, Fixed Arm (18 x 16")	1	1
3		1152604	Kit, Seat Upholstery w/ Hardware, Removable Arm (18 x 16")	1	1
3		1157813	Kit, Seat Upholstery w/ Hardware, Fixed Arm (18 x 16")	1	1
4		1153807	Kit, Arm Pad w/Hardware - Desk Length	1	1
4		1152606	Kit, Arm Pad w/Hardware - Full Length	1	1
5		1152607	Kit, Dust Cap	1	1
6	A	1151514	Kit, Wheel w/ Hardware - Front	1	2
7	B	1080380-NLA	NLA - Wheel Assembly, Composite, Non-Corded Urethane (24")	1	2
7	C,C1	1122726	Wheel Assembly, Composite, Corded Urethane (24")	1	2
8		1063919	Handrim, Composite, 9 PT., Black (24" wheel)	1	2
	D	1065791	Hardware, Handrim Attaching - 9 PT.	1	2
8	G	1055272-NLA	NLA - Handrim, Chrome, 6 PT. (24" Wheel)	1	2
8		1183681	Handrim, Aluminum, 6 PT, Clear Anodized (24" Wheel)	1	2
	D	1042074	Hardware, Handrim Attaching - 6 PT.	1	2
9	C,C1	40117X022	Bearing (1/2" I.D. x 1-3/8" O.D.)	1	2
9	C,C1	40117X024	Bearing (7/16" I.D. x 1-3/8" O.D.)	1	2
10	E	1111924	Wheel Lock Assembly, Push to Lock, Black - Right	1	1
10	E	1111925	Wheel Lock Assembly, Push to Lock, Black - Left	1	1
10	F	1101452	Wheel Lock Assembly, Push to Lock, Black - Right	1	1
10	F	1101453	Wheel Lock Assembly, Push to Lock, Black - Left	1	1
	D	1017677	Kit, Hardware, Wheel Lock (Pair)	1	1
			Bolt, Carriage (1/4-20 x 1-1/4")	1	2
			Nut, Hex Flang (1/4-20)	1	2

Veranda Wheelchair

Footrest packaged with the wheelchair match frame finish. See Section Titled, "Veranda Only" for T93H/T94H Series Front Rigging replacement Options.

At this time, Upholstery replacements for 16" and 20" Chairs are unavailable.

The new front caster used on a 16 and 20" Chair was changed as of March 8, 2011 and for the 18" Chair was changed as of June 21, 2011. "New style caster is available through The Aftermarket Group as TAGRP265093"

NOTE: A - The front caster used on 16" and 20" chairs changed March 8, 2011. The front caster used on the 18" chair changed June 21, 2011. "New style caster is available through The Aftermarket Group as TAGRP265093"
B - No Longer Available. See part #1122726.
C - Prior to 3/8/11, the rear wheel used 1/2" bearings. After 3/7/11, there was a running change in bearing size for the Veranda rear wheel and new wheelchairs could have used 7/16" bearings or 1/2" bearings.
C1 - Part #1122726 will include 7/16" bearings. If you are ordering a wheel replacement, please also order part #40117X022 1/2" bearings (quantity 2 ea. per wheel) so that, if needed, the customer will be able to retrofit the wheel assembly with 1/2" bearings.
D - Not shown.
E - Use for Permanent Arm models. These Wheel Locks are recommended replacements.
F - Use for Removable Arm models. These Wheel Locks are recommended replacements.
G - Chrome Handrim part #1055272 replaced with Clear Anodized Handrim part #1183681.

NOTE: At this time, there are no replacements for the 16 and 20" Back or Seat Upholstery.