



# Footplate Component Kits



## Footplate Component Kits

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1185935	Kit, Footboard, One Piece w/ Hardware (8 - 11" W)	1	1
2			Footboard, One Piece (8-11" W)	1	1
3			Screw (M5 x 16mm)	1	4
4			Nut, Cap (M5)	1	4
1	A	1185936	Kit, Footboard, One Piece w/ Screws (12 - 15" W)	1	1
2			Footboard, One Piece (12-15" W)	1	1
3			Screw (M5 x 16mm)	1	4
4			Nut, Cap (M5)	1	4
5	B	1185937	Kit, Telescoping Tubes w/ Hardware (80 Degree) (Pair)	1	1
6			Tube, Telescoping	1	2
7			Screw, Socket Head Cap (M6 x 1 x 30)	1	2
8			Nut (M6)	1	2
9			Screw (M6 x 1 x 45mm)	1	2
10			Washer (6 x 12 x 1.2)	1	2
11			Plug Button, Black (3/4")	1	2
12			Disk, Locking, Coved	1	4
13	C	1185938	Kit, Telescoping Tubes w/ Hardware (90 Degree) (Pair)	1	1
14			Plug, Footrest Tube	1	2
15			Screw (M6 x 1 x 30mm)	1	2
16			Tube, Support, 90 Degree	1	2
17			Screw, Socket Head (M6 x 50)	1	2
18			Washer (6 x 12 x 1.2)	1	2
19			Washer, Plastic	1	2
20			Cap, Plastic (17 mm)	1	2
21			Disk, Locking, Coved	1	4
22	D	1428812-0601	Kit, Hanger Support w/ Screw	1	1
23			Screw (M8 x 35mm)	1	1
24			Support, Hanger, Elevated Legrest	1	1
25			Indexer, Hanger, 90 Degree - Myon Junior	1	1
26	E	1185932	Kit, Tube, ELR, Upper and Lower (12" W)	1	1
27			Tube, Upper, Elevating Legrest (12" W)	1	1
28			Tube, Lower, Elevating Legrest (12" W)	1	1
29			Screw, CBHC (M8 x 20mm)	1	1
26	E	1185933	Kit, Tube, ELR, Upper and Lower (13 - 15" W)	1	1
27			Tube, Upper, Elevating Legrest (13 - 14" W)	1	1
28			Tube, Lower, Elevating Legrest (13 - 14" W)	1	1
29			Screw, CBHC (M8 x 20mm)	1	1
30	F	1182284	Kit, 3 Point Knob w/ Hardware	1	1
31			Screw, Ball (7.5 mm)	1	1
32			Bushing, Armrest	1	1
33			Nut, Square, Armrest	1	1
34	G	1484580	Kit, Plug Pivot w/ Screw	1	1
35			Pivot, Hanger, 80 Degree	1	1
36			Screw (M6 x 45mm)	1	1
37	H	1185934	Kit, Housing Latch w/ Hardware	1	1

## Footplate Component Kits

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
38			Nut, Lock (M5 x .8)	1	5
39			Bushing, Latch Housing	1	1
40			Housing, Latch	1	1
41			Bracket, Swing, Male	1	1
42			Bracket, Swing, Female	1	1
43			Screw (M5 x .8 x 23mm)	1	3
44			Washer (M6 x 12mm x 0.8mm)	1	1
45			Screw (M5 x .8 x 40mm)	1	1
46			Bushing, Latch	1	1
47			Screw (M5 x .8 x 25mm)	1	1
48			Spring, Double	1	1
48		0572-06603	Spring, Double	1	1
49		0021-60026	Knob, Screw	1	1
50		ZP991451A-NLA	Lever, Footrest Hanger - Right	1	1
50		ZP991450A	Lever, Footrest Hanger - Left	1	1
<p>NOTE: A - Includes items 2, 3, and 4.            B - Includes items 6 - 12.            C - Includes items 14 - 21.            D - Includes items 23, 24, and 25.            E - Includes items 27, 28, and 29.            F - Includes items 31, 32, and 33.            G - Includes items 35 and 36.            H - Includes items 38 - 48.</p>					