


Seat Frame and Hardware

Basic Tilt and Tilt Only or Tilt/Recline and Recline Only



Seat Frame and Hardware

Basic Tilt and Tilt Only or Tilt/Recline and Recline Only

#1 Seat Side Frame and Hardware

Seat Depth	14-15" Wide Right	14-15" Wide Left	16-24" Wide Right	16-24" Wide Left
14"-15" Depth	1072933-NLA	1072934-NLA	N/A	N/A
16"-20" Depth	1072931-NLA	1072932-NLA	1072744-NLA	1072745-NLA
21"-24" Depth	N/A	N/A	1072807-NLA	1072808-NLA

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	Seat Side Frame and Hardware	1	1
2	B	0114702-NLA	NLA - Frame, Seat, Black	1	1
2	C	0208902-NLA	NLA - Frame, Seat, Black	1	1
3		1025736	Locknut (1/4-20)	1	2
4		0050406	Washer, Flat (1/4 x 5/8 x 1/16")	1	4
5		0050506-NLA	NLA - Lockwasher, High Collar, Split (1/4")	1	4
6		0051006-NLA	NLA - Screw, Hex Head Full Thread (1/4-20 x 1-1/2")	1	2
NOTE: A - Includes seat side frame and items 3-6 B - 14" - 20" Wide C - 21" - 24" Wide					