




Armrest Assembly

Non-Flip Back Arms (Before 10/15/98)

Basic Tilt and Tilt Only



14" Deep chairs only.

Armrest Assembly

Non-Flip Back Arms (Before 10/15/98)

Basic Tilt and Tilt Only

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1041522	Arm Tube Assembly w/ Plug Buttons, Full Length, Black - Right/Left	1	1
1		1070859	Arm Tube Assembly w/ Plug Buttons, Desk Length, Black - Right/Left	1	1
2	A	1117328	Kit, Arm Pad w/ Screws, Full Length, Urethane Waterfall, Black (13" L x 2-5/8" W)	1	1
2		1110818	Kit, Arm Pad w/ Screws, Full Length, Urethane Waterfall, Black (13" L x 2-5/8" W) (Pair)	2	1
2	A	1117327	Kit, Arm Pad w/ Screws, Desk Length, Urethane Waterfall, Black (9" L x 2-5/8" W)	1	1
2		1110817	Kit, Arm Pad w/ Screws, Desk Length, Urethane Waterfall, Black (9" L x 2-5/8" W) (Pair)	2	1
3		1027095	Plug Button, Black (7/8")	1	2
4		1055949	Package, Screw, Phillips Pan Head (#10-32 x 1-3/8")	2	1
4	B	1094270-NLA	NLA - Package, Screw, Phillips Pan Head (#10-32 x 1-1/4")	2	1
5		1090768	Kit, Lock Handle w/Rivets (Pair)	2	1
6		47471X074	Spring, Arm Lock	1	1
7		1028594	Pin, Arm Lock, Black	1	1
8		1032245	Grip, Foam Hand (6-7/8" L)	1	1
9		1018225	Spring, Torsion with Button Lock	1	1
10		1028645	Rivet	1	6
11		1064130-NLA	NLA - Guard, Clothing	1	1
12		1035065	Cover, Clothing Guard w/ Pocket, FXDESK & ADJFULL/ADJDESK	1	1
13		0050406	Washer, Flat (1/4 x 5/8 x 1/16")	1	8
14		7930333	Locknut (1/4-20)	25	4
15		1024793-NLA	NLA - Screw, Hex Head (1/4-20 x 1-3/4")	1	4
16		0050506-NLA	NLA - Lockwasher, High Collar, Split (1/4")	1	4
	C	See Chart	Armrest Assembly - 14-15" Deep	1	1

NOTE: A - As of 1/3/02 upholstered waterfall armpads were replaced with black urethane waterfall armpads. Order in pairs with screws for chairs serviced before 1/3/02.
 B - For Black Urethane Waterfall Arm Pads only
 C - Not Shown

Item Number	17	17	18	19	20	21
Description	Arm Tube with Arm Pad Full Length	Arm Tube with Arm Pad Desk Length	Arm Base Assembly	Armrest Socket Front	Armrest Socket Rear	Armrest Mounting Hardware
14" Deep Right	1110821-NLA See Note(s) A	1110819-NLA See Note(s) A	1072943-NLA	1072951-NLA	0121602-NLA	1072949

Armrest Assembly
Non-Flip Back Arms (Before 10/15/98)
Basic Tilt and Tilt Only

Item Number	17	17	18	19	20	21
14" Deep Left	1110822-NLA See Note(s) A	1110820-NLA See Note(s) A	1072944-NLA	1072952-NLA	0121502-NLA	1072950
15" and Deeper Right	1110821-NLA See Note(s) A	1110819-NLA See Note(s) A	1057411-NLA	0121502-NLA	0121302-NLA	1072789-NLA
15" and Deeper Left	1110822-NLA See Note(s) A	1110820-NLA See Note(s) A	1057412	0121602-NLA	0121402-NLA	1072930-NLA

NOTE: A - As of 1/3/02 upholstered waterfall armpads were replaced with black urethane waterfall armpads. Order right and left arm tube with arm pad for chairs serviced before 1/3/02.