


# Headrest Extension Non-Recliner Chairs



## Headrest Extension Non-Recliner Chairs

| Item Number | Note  | Part Number        | Description  | Quantity Per |          |
|-------------|-------|--------------------|--|--------------|----------|
|             |       |                    |  | Package      | Assembly |
| 1           | A,B,F | 8881050181-NL<br>A | NLA - Headrest Extension Assembly, Vinyl, Chrome (14" W)     | 1            | 1        |
| 1           | A,B,F | 8881050182-NL<br>A | NLA - Headrest Extension Assembly, Vinyl, Chrome (15" W)     | 1            | 1        |
| 1           | A,B,F | 8881050183         | NLA - Headrest Extension Assembly, Vinyl, Chrome (16" W)     | 1            | 1        |
| 1           | A,B,F | 8881091224-NL<br>A | NLA - Headrest Extension Assembly, Vinyl, Chrome (16-1/2" W) | 1            | 1        |
| 1           | A,B,F | 8881050184-NL<br>A | NLA - Headrest Extension Assembly, Vinyl, Chrome (17" W)     | 1            | 1        |
| 1           | A,B,F | 8881050185-NL<br>A | NLA - Headrest Extension Assembly, Vinyl, Chrome (18" W)     | 1            | 1        |
| 1           | A,B,F | 8881050186-NL<br>A | NLA - Headrest Extension Assembly, Vinyl, Chrome (19" W)     | 1            | 1        |
| 1           | A,B,F | 8881050187-NL<br>A | NLA - Headrest Extension Assembly, Vinyl, Chrome (20" W)     | 1            | 1        |
| 1           | A,B,F | 8881050188-NL<br>A | NLA - Headrest Extension Assembly, Vinyl, Chrome (22" W)     | 1            | 1        |
| 1           | B,F   | 8881052776-NL<br>A | NLA - Headrest Extension Assembly, Nylon, Chrome (14" W)     | 1            | 1        |
| 1           | B,F   | 8881052777-NL<br>A | NLA - Headrest Extension Assembly, Nylon, Chrome (15" W)     | 1            | 1        |
| 1           | B,F   | 8881052778-NL<br>A | NLA - Headrest Extension Assembly, Nylon, Chrome (16" W)     | 1            | 1        |
| 1           | B,F   | 8881052779-NL<br>A | NLA - Headrest Extension Assembly, Nylon, Chrome (17" W)     | 1            | 1        |
| 1           | B,F   | 8881052780-NL<br>A | NLA - Headrest Extension Assembly, Nylon, Chrome (18" W)     | 1            | 1        |
| 1           | B,F   | 8881052781-NL<br>A | NLA - Headrest Extension Assembly, Nylon, Chrome (19" W)     | 1            | 1        |
| 1           | B,F   | 8881052782-NL<br>A | NLA - Headrest Extension Assembly, Nylon, Chrome (20" W)     | 1            | 1        |
| 2           | B     | 1027095            | Plug Button, Black (7/8")                                    | 1            | 2        |
| 3           | C     | 01058X000-NLA      | NLA - Tube, Headrest Extension                               | 1            | 2        |
| 3           | D     | 1054604-NLA        | NLA - Tube, Headrest Extension                               | 1            | 2        |
| 4           |       | 48811X003          | Tube, Insert   | 1            | 2        |
| 5           |       | 1028753            | Bracket, Headrest, Black - Left                              | 1            | 1        |
| 6           |       | 1028752-NLA        | NLA - Bracket, Headrest, Black - Right                       | 1            | 1        |
| 7           |       | 40163X013-NLA      | NLA - Screw, Phillips Pan Head Tap (#10-32 x 1/2")           | 1            | 2        |
| 8           |       | 1076283            | Screw, Phillips Oval Head w/ Washer (#10-16 x 7/8")          | 1            | 2        |
| 9           | A     | 8881054633-NL<br>A | NLA - Upholstery, Headrest, Vinyl (14" W)                    | 1            | 1        |
| 9           | A     | 8881054634-NL<br>A | NLA - Upholstery, Headrest, Vinyl (15" W)                    | 1            | 1        |
| 9           | A     | 8881054638-NL<br>A | NLA - Upholstery, Headrest, Vinyl (16" W)                    | 1            | 1        |
| 9           | A     | 8881054635-NL<br>A | NLA - Upholstery, Headrest, Vinyl (17" W)                    | 1            | 1        |
| 9           | A     | 8881054637-NL<br>A | NLA - Upholstery, Headrest, Vinyl (18" W)                    | 1            | 1        |

## Headrest Extension Non-Recliner Chairs

| Item Number  | Note | Part Number        | Description                                   | Quantity Per |          |
|--|------|--------------------|---|--------------|----------|
|  |      |                    |   | Package      | Assembly |
| 9  | A    | 8881054636-NL<br>A | NLA - Upholstery, Headrest, Vinyl (19" W)     | 1            | 1        |
| 9  | A    | 8881054639-NL<br>A | NLA - Upholstery, Headrest, Vinyl (20" W)     | 1            | 1        |
| 9  | A    | 8881054640-NL<br>A | NLA - Upholstery, Headrest, Vinyl (22" W)     | 1            | 1        |
| 9  |      | 8881052586-NL<br>A | NLA - Upholstery, Headrest, Nylon (14" W)     | 1            | 1        |
| 9  |      | 8881052587-NL<br>A | NLA - Upholstery, Headrest, Nylon (15" W)     | 1            | 1        |
| 9  |      | 8881052588-NL<br>A | NLA - Upholstery, Headrest, Nylon (16" W)     | 1            | 1        |
| 9  |      | 8881091223-NL<br>A | NLA - Upholstery, Headrest, Vinyl (16-1/2" W) | 1            | 1        |
| 9  |      | 8881052589-NL<br>A | NLA - Upholstery, Headrest, Nylon (17" W)     | 1            | 1        |
| 9  |      | 8881052590-NL<br>A | NLA - Upholstery, Headrest, Nylon (18" W)     | 1            | 1        |
| 9  |      | 8881052591-NL<br>A | NLA - Upholstery, Headrest, Nylon (19" W)     | 1            | 1        |
| 9  |      | 8881052592-NL<br>A | NLA - Upholstery, Headrest, Nylon (20" W)     | 1            | 1        |
|  | E,G  | 1098351-NLA        | NLA - Instruction Sheet                       | 1            | 1        |
|  | C,E  | 47447X009-NLA      | NLA - Screw, Phillips (#8 x 3/8")             | 1            | 6        |
| <p>NOTE: A - Specify upholstery color when ordering, see category titled, "Upholstery Color Options."<br/>           B - Includes items 2-9.<br/>           C - For use on a Headrest Extension Tube with snap fasteners (Before 6/95).<br/>           D - For use on a Headrest Extension Tube without snap fasteners (After 5/95).<br/>           E - Not Shown.<br/>           F - Headrest extension assembly not available in 71P Silver Vein color frame finish.<br/>           G - Item is no longer available from the Invacare Service Parts Department. Refer to the Invacare Website under "Product Literature Search". Enter the part number in the "Keywords" field and download a free PDF version copy.</p> |      |                    |   |              |          |