



# Frames without Hardware



# Frames without Hardware

## Standard

Side	Chrome	Color
10" Deep - Right	1100935-NLA	8881100933
	See Note(s) B1, B	See Note(s) B, B1, A
10" Deep - Left	1100936-NLA	8881100934
	See Note(s) B1, B	See Note(s) B1, B, A
12" Deep - Right	1100939-NLA	8881100937
	See Note(s) B1, B	See Note(s) B1, B, A
12" Deep - Left	1100940-NLA	8881100938
	See Note(s) B1, B	See Note(s) B, A, B1
14" Deep - Right	1100943-NLA	8881100941
	See Note(s) B1, B	See Note(s) B, B1, A
14" Deep - Left	1100944-NLA	8881100942
	See Note(s) B, B1	See Note(s) B, B1, A

NOTE: A - Specify frame finish, see category titled "Frame Finishes, Packaging, Literature and Upholstery Colors".  
 B - Hinge Pin Replacement Hardware Kit located on the Frame Hardware page is used on Frames w/ welded-on Hinge Pins manufactured before 4/12/13. Side Frames after 4/11/13 have Hinge Pin Saddle Plates w/ square holes to attach REMOVABLE Non-Welded Hinge Pins.  
 B1 - All Side Frames ordered after 4/11/13 require new Side Frame Hardware that will include the removable Hinge Pins also. Removable Hinge Pins are not sold separately.

## Reclining

Side	Chrome	Color
10" Deep - Right	1100947-NLA	8881100945
	See Note(s) B, B1	See Note(s) A, B, B1
10" Deep - Left	1100948-NLA	8881100946
	See Note(s) B1, B	See Note(s) A, B, B1
12" Deep - Right	1100951-NLA	8881100949
	See Note(s) B1, B	See Note(s) B1, B, A
12" Deep - Left	1100952-NLA	8881100950
	See Note(s) B1, B	See Note(s) B, A, B1
14" Deep - Right	1100955-NLA	8881100953
	See Note(s) B1, B	See Note(s) B1, B, A
14" Deep - Left	1100956-NLA	8881100954
	See Note(s) B1, B	See Note(s) B1, B, A

NOTE: A - Specify frame finish, see category titled "Frame Finishes, Packaging, Literature and Upholstery Colors".  
 B - Hinge Pin Replacement Hardware Kit located on the Frame Hardware page is used on Frames w/ welded-on Hinge Pins manufactured before 4/12/13. Side Frames after 4/11/13 have Hinge Pin Saddle Plates w/ square holes to attach REMOVABLE Non-Welded Hinge Pins.  
 B1 - All Side Frames ordered after 4/11/13 require new Side Frame Hardware that will include the removable Hinge Pins also. Removable Hinge Pins are not sold separately.

# Frames without Hardware

## Adjustable Angle Back

Side	Chrome	Color
10" Deep - Right	1100959-NLA	8881100957
	See Note(s) B1, B	See Note(s) B1, B, A
10" Deep - Left	1100960-NLA	8881100958
	See Note(s) B1, B	See Note(s) B1, B, A
12" Deep - Right	1100963-NLA	8881100961
	See Note(s) B1, B	See Note(s) B1, B, A
12" Deep - Left	1100964-NLA	8881100962
	See Note(s) B1, B	See Note(s) B1, B, A
14" Deep - Right	1100967-NLA	8881100965
	See Note(s) B1, B	See Note(s) B1, B, A
14" Deep - Left	1100968-NLA	8881100966
	See Note(s) B1, B	See Note(s) B1, A, B

NOTE: A - Specify frame finish, see category titled "Frame Finishes, Packaging, Literature and Upholstery Colors".  
 B - Hinge Pin Replacement Hardware Kit located on the Frame Hardware page is used on Frames w/ welded-on Hinge Pins manufactured before 4/12/13. Side Frames after 4/11/13 have Hinge Pin Saddle Plates w/ square holes to attach REMOVABLE Non-Welded Hinge Pins.  
 B1 - All Side Frames ordered after 4/11/13 require new Side Frame Hardware that will include the removable Hinge Pins also. Removable Hinge Pins are not sold separately.