


Back Canes - SX5



Back Canes - SX5 Only

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1110018-NLA	NLA - Back Cane Assembly, Chrome - Right/Left (16" H)	1	1
1	A,B	8881117026	Back Cane Assembly, Color - Right/Left (16" H)	1	1
1	A	1110020-NLA	NLA - Back Cane Assembly, Chrome - Right/Left (18" H)	1	1
1	A,B	8881181863	Back Cane Assembly, Color - Right/Left (18" H)	1	1
2		1021261	Grip, Hand (6")	1	1
3	C	1110008-NLA	NLA - Cane, Back, Chrome - Right/Left (16" H)	1	1
3	C	1110010-NLA	NLA - Cane, Back, Chrome - Right/Left (18" H)	1	1
4		1025195	Locknut (1/4-20)	1	1
5		1056551	Insert	1	1
6		1057955	Screw, Hex Head (1/4-20 x 1-1/4")	1	1

NOTE: A - Includes items 2-6.
 B - Specify 71P Silver Vein Frame finish when ordering.
 C - Back Cane is not available in color, order item #1.