



Clamp-On Style Wheel Locks

Removable Arm, SX5RC w/ 20" or 22" Rear Wheels only



Clamp-On Style Wheel Locks

Removable Arm, SX5RC w/ 20" or 22" Rear Wheels only

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	E,F	1080332-NLA	NLA - Wheel Lock Assembly, Black - Right	1	1
1	E,F	1080333-NLA	NLA - Wheel Lock Assembly, Black - Left	1	1
1	D,E	8881181914	Wheel Lock Assembly, Color - Right	1	1
1	D,E	8881181915	Wheel Lock Assembly, Color - Left	1	1
2		47376X024	Tip, Rubber Wheel Lock Handle, Black	1	1
3	A	1030987	Hardware, Wheel Lock Attaching	1	1
4			Bolt, Lug (5/16-24 x 1-3/8")	1	2
5			Locknut (5/16-24)	1	2
	B,C	1048656	Kit, Wheel Lock Shoe Spacer	1	1
			Spacer, Wheel Lock Shoe (1-45/64")	1	1
			Screw, Set Socket Cup Point (#8-32 x 3/16")	1	2

NOTE: A - Includes items 4 and 5.
 B - Kit quantity covers one side only and used with 20" rear wheels only.
 C - Not Shown
 D - Specify 71P Silver Vein Frame Finish when ordering.
 E - Includes item 2.
 F - Chrome is No Longer Available. Order in quantity of 2 to match