


Frame Deck Assembly



Frame Deck Assembly

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1165043	Kit, Head Deck w/ Hardware	1	1
2		1165046	Kit, Thigh and Foot Deck w/ Hardware	1	1
3		1149706-NLA	NLA - Frame, Bed	1	1
4		1147550	Ratchet, Knee	1	2
5		1147565	Clevis Pin, Deck to Frame (36" W)	1	4
6		PIN51634-NLA	NLA - Pin, Clevis (5/16" x 3/4")	1	2
7		1186877-NLA	NLA - Bolt, Attachment (1/4-20)	1	2
8		1186909-NLA	NLA - Nut, Push (5/16)	1	2
9		1186857-NLA	NLA - Washer, Nylon (3/8 x 3/4 x 1/16")	1	2
10		1186860-NLA	NLA - Lock Nut, Purple Insert (1/4-20)	1	2
11		1159596	Kit, Bracket, Bed End - Foot (Pair)	1	1
			Bracket, Bed End, Adjustable, Foot - Right	1	1
			Bracket, Bed End, Adjustable, Foot - Left	1	1
12		1161887	Kit, Bed End Bracket, Head (Pair)	1	1
			Bracket, Bed End, Adjustable, Head - Right	1	1
			Bracket, Bed End, Adjustable, Head - Left	1	1