



Rear Wheel, Handrim and Mounting Hardware



SSM_020A

Rear Wheel, Handrim and Mounting Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1152229	Kit, Wheel 24", 7/16" Axle, Composite Urethane - 9 PT. Plastic Coated Handrim	1	1
2			Wheel Assembly, Composite, Corded Urethane (24")	1	1
3			Handrim, Composite, 9 PT., Black (24" wheel)	1	1
4			Hardware, Handrim Attaching - 9 PT.	9	1
9	B	1570524	Kit, Rear Wheel Mounting Hardware	1	1
10			Screw, Hex Head (7/16-20 x 4-1/2")	1	1
11			Locknut (7/16-20)	1	1
NOTE: A - Includes items 2-4. B - Includes items 10 and 11.					