


Drive Wheels, Casters and Forks



SCP_007A

Drive Wheels, Casters and Forks

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1186182	Kit, Drive Wheel FF Assembly w/ Hardware Includes: (1) ea. Drive Wheel (1) ea. M10 Nut (1) ea. Washer	1	1
2		1186178	Kit, Fork, Caster w/ Hardware Includes: (1) ea. Fork (1) ea. Fork Stud (1) ea. Nylon Nut	1	1
3		1186179	Kit, Solid Caster Wheel Assembly Includes: (1) ea. Wheel Assembly (2) ea. Washer (1) ea. Screw (1) ea. Nylon Nut	1	1