


# Seat and Seat Base



SCP\_020A

## Seat and Seat Base

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1186188	Seat Assembly Includes: (1) ea. Rear Base Plate (1) ea. "L" Type Release Lever (1) ea. Right Lever Tip, Black (1) ea. Spring Lever (3) ea. Knob (1) ea. Square Plug (1) ea. Right Locating Lever (1) ea. Set Screw (1) ea. Spacer (1) ea. Hex Screw (1) ea. Spring Washer (1) ea. Push Tie-mount (1) ea. Headrest Assembly (1) ea. Headrest Cover, Black (1) ea. Backrest Cover, Black (1) ea. Seat Cover, Black	1	1
2	A	1186189	Kit, Replacement Knob Includes: (3) ea. Knob	1	1
		1186351	Kit, Seat Refurbish Includes: 2 (ea.) Armrest Pad 2 (ea.) Knob 4 (ea.) Countersunk Screws 1 (ea.) Headrest Cover 1 (ea.) Backrest Cover 1 (ea.) Seat Cover 1 (ea.) Instruction Sheet	1	1
NOTE: A - Not shown on illustration.					