


Footplate and Footplate Bracket



SCP_015A

Footplate and Footplate Bracket

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1186194	Kit, Footplate w/ Hardware Includes: (1) ea. Footplate (1) ea. Brass Sleeve (1) ea. Hex Head Bolt (2) ea. Washer (1) ea. Nylon Nut	1	1
2		1186195	Kit, Footplate Bracket w/ Hardware Includes: (1) ea. Screw (1) ea. Nut (1) ea. Footplate Bracket (2) ea. Screw (2) ea. Nylon Nut	1	1