


A.B.S. Curved Back Assembly

Basic Recline (4XP)



Material Back Width x Height	Overall Back Height	Item #1 Back Pan	Item #2 1" Back Foam	Item #3 Back Cover	Item #4 Back Asbly w/Fixed Hdwr
16" x 20"	24" Fixed	0125319-NLA	0009820-NLA	0011321-NLA	1072873-NLA
18" x 20"	24" Fixed	0125519-NLA	0010220-NLA	0011721-NLA	1072875-NLA

NOTE: Some foam size part numbers can accommodate two width and two back heights by turning the foam around 90 degrees. Example would be a (16" W x 20" H) can be (20" W x 16" H) in other applications.

A.B.S. Curved Back Assembly

Basic Recline (4XP)

Basic Recline (4XP)-MSR

Material Back Width x Height	Overall Back Height	Item #1 Back Pan	Item #2 1" Back Foam	Item #3 Back Cover	Item #4 Back Asbly w/Sliding Hdwr
16" x 14"	21"-22"H	0129019-NLA	1073070	0124821-NLA	1077097
18" x 14"	21"-22"H	0129219-NLA	1073071-NLA	0125021-NLA	1077099-NLA
16" x 16"	23"-25"H	0049019-NLA	0009720	0011221-NLA	1072878
18" x 16"	23"-25"H	0049219-NLA	0010120-NLA	0011621-NLA	1072880-NLA

NOTE: Some foam size part numbers can accommodate two width and two back heights by turning the foam around 90 degrees. Example would be a (16" W x 20" H) can be (20" W x 16" H) in other applications.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1-4		See Chart	#s 1-4 Basic Recline (4XP)	1	1
1-4		See Chart	#s 1-4 Basic Recline (4XP)-MSR	1	1
5		0029002-NLA	NLA - Bracket, Non-Sliding Back Pan Mounting	1	4
6		7930323-NLA	NLA - Package, Locknut (1/4-20)	25	4
7		1029545	Screw, Socket Flat Head (1/4-20 x 3/4")	1	4
8		0047800-NLA	NLA - Block, Glide	1	4
9		0052906-NLA	NLA - Screw, Socket Button Head (1/4-20 x 3/4")	1	4
	A,D	1125723	Kit, Hook Fastener w/ Adhesive Back	1	1
	A,B	0230200-NLA	NLA - Flap, Back Pan - Tarsys	1	1
	A,C	0230300-NLA	NLA - Flap, Back Pan - Tarsys	1	1

NOTE: A - Not Shown
B - 16" Wide
C - 18" Wide
D - Hook Fastener Adhesive Kit is used to attach the Back Cover to the Back Pan and must be ordered when a new Back Pan or Back Assembly is requested.