


Forks and Hardware



Forks and Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1033202	Kit, Fork Hardware (Pair)	1	1
2			Cap, Head Tube, Black	1	2
3			Locknut (1/2-20)	1	2
4			Washer, Nylon (5/8 x 1 x 3/32")	1	2
9		1001801	Bearing, Flange (1/2")	1	2
10	B	1056640	Fork and Stem Assembly, Black (8" & 6 x 2" Caster)	1	1
11	C	1111305	Fork and Stem Assembly, Steel, Black (8")	1	1
NOTE: A - Includes Items 2 thru 4. B - Used for the Veranda (K0004) model C - Used for the Veranda (K0003) model					