


# Midline Assembly



## Midline Assembly

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		0063100-NLA	NLA - Collar, 2 Holes	1	1
2		0223506-NLA	NLA - Screw, Socket, Cap, Button Head (1/4-20 x 1/2")	1	1
3		1077307-NLA	NLA - Screw, Set (1/4-20 x 1/4")	1	4
4		0063000-NLA	NLA - Spring, Assist	1	1
5		1077306-NLA	NLA - Bearing, Nylon (1/2")	1	2
6		0063402-NLA	NLA - Plate, Mounting	1	1
7		0191001-NLA	NLA - Bar, Mounting	1	1
8		7930304-NLA	NLA - Screw, Socket Head (1/4-20 x 1-1/4")	25	3
9		0063202-NLA	NLA - Weldment, Collar	1	1
10		1077308-NLA	NLA - Screw, Set (1/4-20 x 1/2")	1	1
11		0053200	Rivnut (1/4-20 x 11/16")	1	1
12		1077321-NLA	NLA - Knob, Tee (1/4-20 x 1-3/4")	1	1
13		0267606-NLA	NLA - Spacer (17/64 x 21/64 x 1/2")	1	1
14		0050100-NLA	NLA - Washer, Nylon Flat (1/4 x 1/2 x 1/16")	1	2
15		0064600-NLA	NLA - Spring, Compression	1	1
16		7930333	Locknut (1/4-20)	25	1
17		1077309-NLA	NLA - Screw, Set (1/4-20 x 1")	1	1
18		0312524-NLA	NLA - Plate, Mounting, MKIV Std/ PSF/ PSR Joystick	1	1
18		1077310-NLA	NLA - Plate, Mounting, MCC MKIII Std. Joystick	1	1
18		1077311-NLA	NLA - Plate, Mounting, 1556 H.D. / 1810M Short Std. Joystick	1	1
18		1077312-NLA	NLA - Plate, Mounting MCCX MKIII Joysticks	1	1
18		0314602-NLA	NLA - Plate, Mounting, P.S. Joysticks	1	1
19		1077249-NLA	NLA - Screw, Socket Head (#10-32 x 3/4")	1	2
20		0068000-NLA	NLA - Screw, Phillips Flat Head (#6 x 5/8")	1	4
21		0106600-NLA	NLA - Screw Phillips Flat Head (#8-32 x 1/2")	1	2
22		1025716	Screw, Phillips Pan Head (#10-32 x 11/16")	1	2
23		0224806-NLA	NLA - Washer (#10)	1	2
24		1077313-NLA	NLA - Mid-Line Assembly, Less Half Clamps, MKIV Std./ PSF/ PSR Joystick	1	1
24		1077314-NLA	NLA - Mid-Line Assembly, Less Half Clamps, MCC MKIII Std. Joystick	1	1
24		1077315-NLA	NLA - Mid-Line Assembly, Less Half Clamps, 1556 Heavy Duty/ 1810M Short Std. Joystick	1	1
24		1077316	Mid-Line Assembly, Less Half Clamps, MCCX MKIII Joystick	1	1
24		1077317-NLA	NLA - Mid-Line Assembly, Less Half Clamps, P.S. Joystick	1	1
25		1027266	Half Clamp, Non-Threaded Holes	1	1
26		1027267-NLA	NLA - Half Clamp, Threaded Holes	1	1