


# Mechanical Legrest Actuator Linkage Without M.S.R. Basic Recline (4XP)



## Mechanical Legrest Actuator Linkage Without M.S.R. Basic Recline (4XP)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1024793-NLA	NLA - Screw, Hex Head (1/4-20 x 1-3/4")	1	1
2		0050100-NLA	NLA - Washer, Nylon Flat (1/4 x 1/2 x 1/16")	1	7
3		0233222-NLA	NLA - Spacer, Plastic	1	1
4		0050406	Washer, Flat (1/4 x 5/8 x 1/16")	1	2
5		7930333	Locknut (1/4-20)	25	3
6		0055506-NLA	NLA - Screw, Button Head, Socket (5/16-18 x 2")	1	2
7		0060300-NLA	NLA - Bushing, Bronze (5/16 x 7/16 x 1-1/4")	1	1
8		0055800-NLA	NLA - Washer, Nylon (7/16 x 3/4 x 3/32")	1	2
9		0054906	Locknut, High Profile (5/16-18)	1	2
10		0051606-NLA	NLA - Screw, Cap, Button Head, Socket (1/4-20 x 1")	1	1
11		0217324-NLA	NLA - Link, Small, Black	1	1
12		0233506-NLA	NLA - Screw, Hex Head (1/4-20 x 2-1/2")	1	1
13		0217024-NLA	NLA - Bushing, Back Post/Leg Link	1	1
14		0217124-NLA	NLA - Link, Back Post Leg, Black	1	1
15		0220924-NLA	NLA - Spacer, Rod End	1	1
16		0083800-NLA	NLA - Rod, End, Male (5/16-24)	1	1
17		0141405-NLA	NLA - Nut, Hex (5/16-24)	1	1
18	B	0035902-NLA	NLA - Rod, Push, Rear Tube, Short, Black (12-1/2")	1	1
18	C	0222502-NLA	NLA - Rod, Push, Rear Tube, Long, Black (14")	1	1
19	D	0035702-NLA	NLA - Rod, Push, Front Tube, AT5544 Legrest, Black	1	1
19	E	0036402	NLA - Rod, Push, Front Tube, Genius Legrest, Black	1	1
20		0026400-NLA	NLA - Receptacle, Power Legrest	1	1
21		0053306-NLA	NLA - Screw, Allen Head (1/4-28 x 1")	1	1
22	A	1077510-NLA	NLA - Mechanical Legrest Hardware Assembly	1	2

NOTE: A - Less items 18, 19  
B - Short 12-1/2"  
C - Long 14"  
D - AT5544 Legrest  
E - Genius Legrest