


Anti-Tip Wheel Extension Assembly

Vent Tray Option Only (After 6/30/98)

2nd Generation Suspension Arms



Anti-Tip Wheel Extension Assembly

Vent Tray Option Only (After 6/30/98)

2nd Generation Suspension Arms

Anti-Tip for Rubber Element

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1070397-NLA	NLA - Screw, Hex Head (5/16-18 x 3-1/4")	1	2
2		1079782-NLA	NLA - Spacer (21/64 x 1/2 x 3/8")	1	4
3		1026158	Locknut (5/16-18)	1	10
4		1065215	Screw, Socket Button Head w/ Patch (5/16-18 x 3/4")	1	8
5		1079778-NLA	NLA - Bracket, Extension, Black	1	4
6		1032190-NLA	NLA - Wheel Assembly (3" x 1")	1	2
7		1039860	Screw, Hex Head w/ Patch (5/16-18 x 3/4")	1	6
8		1028769	Locknut (5/16-18)	1	6
9		1079780-NLA	NLA - Anti-Tipper Extension Assembly	1	1

Anti-Tip for Shock Assembly

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
10		1070397-NLA	NLA - Screw, Hex Head (5/16-18 x 3-1/4")	1	4
11		1079782-NLA	NLA - Spacer (21/64 x 1/2 x 3/8")	1	4
12		1026158	Locknut (5/16-18)	1	8
13		1065215	Screw, Socket Button Head w/ Patch (5/16-18 x 3/4")	1	4
14		1079778-NLA	NLA - Bracket, Extension, Black	1	4
15		1032190-NLA	NLA - Wheel Assembly (3" x 1")	1	2
16		1079781-NLA	NLA - Anti-Tipper Extension Assembly	1	1