



# Tub Bars


**Model 100-3 - Chrome**  
**Model 110-3 - Vinyl (Shown)**  
**Model 1816 - Narrow Wall, Chrome**



**Model 725-3 - Vinyl**




**Model 1200-3 - Chrome**  
**Model 1210-3 - Vinyl**



## Tub Bars

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1073204-NLA	NLA - Wrench, Mounting Hex Nut (1/2")	1	1
2		1073208-NLA	NLA - Bolt, Hex Cap (1/4-20 x 3")	1	2
3		1073517	Bolt, Hex (1/4-20 x 1-3/4")	1	2
4		16720X001-NLA	NLA - Bolt, Hex Cap (1/4-20 x 2")	1	2
5		1072003-NLA	NLA - Pad, Non-Skid Adhesive Back	1	4
6		1072004	Pad, Non-Skid Adhesive Back	1	4

# Tub Bars



## Tub Bars

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1052984-NLA	NLA - Knob with Threaded Rod, White	1	1
2		1052959-NLA	NLA - Pad, Non-Skid without Holes	1	1
3		1052960-NLA	NLA - Pad, Non-Skid without Holes	1	1